

Patient Monitor Mounting

Quick Installation and Maintenance Instructions



PM5-

Medical Facility Responsibility

Preventive maintenance checks must be performed regularly to maintain the quality and performance of this product. Any parts that may be broken, missing, worn, distorted, or contaminated in any way should not be used and all affected parts should be replaced immediately. Should the necessity of any repair be suspected; please contact your local Amico distributor.



WARNING: It is the responsibility of the end user to ensure all aspects of the installation are covered by following the full manual found at: <https://www.amico.com/sites/default/files/product/downloads/amico-aa-5in-slide-in-shelf-with-vesa-plate-manual.pdf>.



NOTE: Only the *MONITOR* rail channel (-MON) set up is shown. For other rail type configurations, please refer to the full manual.

Installation Tools



Hex Key: 1/8"

Installation to Mounting Platform (MRS) - Installation to Monitor Rail Channel



NOTE: The *MONITOR BRACKET STOPPER (TOP)* and the *MONITOR BRACKET STOPPER (BOTTOM)* are there to prevent the *MOUNTING ADAPTER* from sliding out of the *MRS RAIL*. The *MOUNTING ADAPTER* must be taken out of the top of the *MRS RAIL* with the *MONITOR BRACKET STOPPER (TOP)* pressed back into the *MRS RAIL*. (Figure 1)

1. Lift up the *SHELF* and guide the *MOUNTING ADAPTER* into the top of *MRS RAIL*. (Figure 1)
2. Once the *SHELF* is in the desired position on the *MRS RAIL*, tighten 2x *BLACK SET SCREW CUPS* to the tightest possible position on the *MOUNTING ADAPTER* using a 1/8" Hex key.






Figure 1

Installation to Device - Installation to Patient Monitor with VESA Pattern

1. Place *MONITOR* face down on a padded surface and attach *VESA PLATE* to *MONITOR* using 4x *M4 SCREWS*. (Figure 2)
2. Slide *VESA PLATE* into *SHELF* (Figure 3). Note that the *SHELF* should be securely attached to the *MOUNTING PLATFORM* at this stage and the *PLUNGER PIN* should be disengaged. To disengage the *PLUNGER PIN*, see instructions in "Mounting to Platform" section found in the full manual.
3. Pull *PLUNGER PIN* slightly back. Rotate 90° again and release. This should re-engage the *PLUNGER PIN* and allow protrusion through the *SHELF* and into the bottom slot on the *VESA PLATE*. If the *PLUNGER PIN* does not engage, pull back, rotate slightly and retry. Figure 3 shows an exploded version (NOTE: the *MRS RAIL* has been omitted for clarity).
4. Tighten 4x *PLASTIC SCREWS* on the *VESA PLATE*. (Figure 4)



Patient Monitor Mounting

Quick Installation and Maintenance Instructions




PM5-

Adjustments

The **MONITOR** can be rotated to a landscape/portrait orientation. To achieve this, follow the steps below.

1. Fully loosen 4x **PLASTIC SCREWS**, disengage the **PLUNGER PIN** by pulling it back and rotating 90°. (**Figure 5**)
2. Carefully slide the **MONITOR** (attached to **VESA PLATE**) from the **SHELF**. (**Figure 6**)
3. Rotate the **MONITOR** to the desired orientation and slide the **VESA PLATE** back into the **SHELF**.
4. Follow steps 3 and 4 from the Mounting Devices section in the full manual to complete adjustment.



Maintenance

It is recommended that once every three (3) months, all locks, bolts, and screws are visually inspected for signs of product wear and tear or damage. Check the product for looseness in the mount or orientation, and ensure the adapter is secured to the rail or wall channel system. Ensure the (2) set screws on the channel and (2) thumb screws are fastened to the tightest possible position and ensure device is secure.

Warranty

During the term of warranty: Within the first twelve (12) months from the date of shipment, Amico Accessories will repair or replace any part which is proven to be defective at no cost. After the twelve (12) month period, Amico Accessories will send the parts to the customer free of charge, however shipping and installation will be borne by the customer.

The warranty is valid only when the product has been properly installed according to Amico Accessories specifications, used in a normal manner, and serviced according to factory recommendations. It does not cover failures due to damage which occur in shipments or failures which resulted from accidents, misuse, abuse, neglect, mishandling, alteration, misapplication or damage that may be attributable to Force Majeure.

AMICO ACCESSORIES DOES NOT HONOR VERBAL STATEMENTS CONCERNING THE WARRANTY.

The distributor and/or dealer are not sanctioned to create verbal warranties about the product described in this agreement. Any statements will not be honored or be made part of the agreement of sale. This document is the final complete and exclusive terms of the agreement.

THIS WARRANTY IS INCLUSIVE AND REPLACES ALL OTHER WARRANTIES.

Amico Accessories shall not, under any circumstances be liable for incidental or consequential damages including, but not limited to, profit, loss of sales or injuries to person(s) or property.

Correction of non-compliance as noted above will result in completion of all liabilities of Amico Accessories whether based on agreement, neglect or changed materials, designs or specifications without notice.

All claims for warranty must first be approved by Amico Accessories Customer Service Department: (info@amico-accessories.com or 1-877-264-2697). A valid Return Goods Authorization number must be obtained from Amico Accessories prior to commencement of any warranty claim.