

# Roll Stand

## Quick Installation and Maintenance Instructions



RS\_-42\_-PKG/RS\_-LCD-PKG

### Medical Facility Responsibility

Preventive maintenance checks must be performed regularly to maintain the quality and performance of this product. Any parts that may be broken, missing, worn, distorted or contaminated in any way should not be used and all affected parts should be replaced immediately. Should the necessity of any repair be suspected; please contact your local Amico distributor.



**WARNING:** It is the responsibility of the end user to ensure all aspects of the installation are covered by following the full manual found at:

<https://www.amico.com/sites/default/files/product/downloads/amico-aa-rst-roll-stand-manual.pdf>

### Installation Tools




3/16" Hex Key



#2, #3 Phillips Head Screwdrivers



1/2" (12.7 mm) Socket Wrench





**WARNING:** If the product was purchased with counterweight, perform Installing Counterweight step 1 and/or step 2 from the full manual. If not, proceed to the Installing Casters & Post into Base section.

### Installing Casters & Post into Base




**WARNING:** Do not use a hammer on the *CASTERS* or *ROLL STAND* base when installing *CASTERS*. Doing so will damage the *CASTERS*.

1. Push each of the 5x *CASTERS* into a leg by hand until it clicks in place, as shown in **Figure 1**.
2. Insert the *POST* into the *BASE ADAPTER*. Using a 1/2" socket wrench, secure by inserting the *SCREW* through the *LARGER WASHER*, *SMALLER WASHER*, the *STAND BASE* underside, and *BASE ADAPTER* into the *POST*, as seen in **Figure 2**.



### Locking & Unlocking Casters

1. To lock the *CASTERS*, step on the *LOCKING LEVER* such that it clicks in the down position.
2. To unlock the *CASTERS*, lift the *LOCKING LEVER* up as seen in **Figure 3**.



# Roll Stand

## Quick Installation and Maintenance Instructions



RS\_-42\_-PKG/RS\_-LCD-PKG

### Locking & Unlocking Casters

1. Using a 1/2" socket wrench, insert the **SCREW** through the **BRONZE WASHER**, the **NYLON WASHER**, the **MONITOR HEAD**, the black **PLASTIC BUSHING**, and the **BASE ADAPTER** into the **POST** as seen in **Figure 4**.
2. Cover the **HEAD** of the **SCREW** with the Amico doming label by placing the adhesive side down onto the **MONITOR HEAD** over the **SCREW**.




Figure 4

### Troubleshooting Guide

| Symptom                             | Possible Cause                                     | Solution                                    |
|-------------------------------------|--|---|
| Roll Stand is difficult to maneuver | Locked or unsecured caster                         | Refer to instructions (Adjustments)         |
| Device is loose on Roll Stand head  | Device secured improperly/lock plunger not engaged | Refer to instructions (Device Installation) |

### Maintenance

It is recommended that once every three (3) months, all locks, bolts and screws are visually inspected for signs of product wear and tear or damage.

### Warranty

Amico Accessories Inc. warrants all mounting accessories to be free from defects in material and workmanship for a period of twelve (12) months from the date of shipment. Within this period Amico Accessories Inc. will repair or replace any part which is proven to be defective.

Amico Accessories Inc. will warrant its materials to be free from defect for an additional period of four (4) years, (five [5] years from the date of shipment). Within this period, Amico Accessories Inc. will replace any part which is proven to be defective, at no charge. Shipping and Installation costs after the first twelve (12) months will be borne by the Customer.

This warranty is valid only when the product has been properly installed according to Amico Accessories Inc. specifications, used in a normal manner and serviced according to factory recommendations. It does not cover products that are not manufactured by Amico Accessories Inc. It does not cover failures due to damage which occurs in shipments or failures which result from accidents, misuse, abuse, neglect, mishandling, alteration, misapplication or damage that may be attributable to acts of God.

Amico Accessories Inc. shall not be liable for incidental or consequential damages resulting from use of equipment.

All claims for warranty must first be approved by Amico Accessories Inc. A valid Return Goods Authorization (RGA) number must be obtained from Amico Accessories Inc. prior to commencement of any service work. Warranty work, which has not been pre-authorized by Amico Accessories Inc., will not be reimbursed.

#### **AMICO ACCESSORIES INC. DOES NOT HONOR VERBAL STATEMENTS CONCERNING THE WARRANTY.**

The distributor and/or dealer are not sanctioned to create verbal warranties about the product described in this agreement. Any statements will not be honored or be made part of the agreement of sale. This document is the final complete and exclusive terms of the agreement.

#### **THIS WARRANTY IS INCLUSIVE AND REPLACES ALL OTHER WARRANTIES.**

Amico Accessories Inc. shall not, under any circumstances be liable for incidental or consequential damages including, but not limited to, profit, loss of sales or injuries to person(s) or property. Correction of noncompliance as noted above will result in completion of all liabilities of Amico Accessories Inc., whether based on agreement, neglect or changed materials, designs or specifications without notice.

All claims for warranty must first be approved by Amico Accessories Inc. Service Department: accessories@amico.com or 1.877.264.2697. A valid Return Goods Authorization number must be obtained from Amico Accessories Inc. prior to commencement of any warranty claim.