Retro-Fit Manual Diamond 3 Outlet





Important Note

Due to unknown factors of the existing equipment, projects requiring Retro-Fit Conversions must be fully investigated by the installer.

It is the responsibility of the installer to determine if the Retro-Fit will function correctly and if additional components are required. Amico is not responsible for any additional components required.

Retro-Fit Instructions

This Retro-Fit should only take about 10 minutes to install. There is no need to shut off the gas. Please note that this application only Retro-Fits the front assembly and no parts will be sent along for the back rough-in. If you require a Retro-Fit to change the back rough-in you'll need to order the Diamond 2 Retro-Fit.

Ohmeda Diamond 2 & 3 Outlet Remove screws from front Facia.



Amico has made this latch valve compatible to the existing rough-in. The barrel is shorter and the diameter is wider. Note the indexing of the pin is on the bottom.



Front View



Side View

Retro-Fit Instructions

Screw in the Amico Diamond 3 Latch Valve and 4.5" chrome frame.



A completed Diamond 3 assembly conversion from Ohmeda to Amico.

MEDICAL AIR MEDICAL AIR MEDICAL AIR MICO

WARNING

Ohmeda issued a Mandatory Field modification for its back body on Diamond 3 Outlets manufactured between 1990 and 1993. If this Ohmeda field modification was not applied, the Outlet MUST NOT be Retrofitted since this could cause total loss of gas flow through the Outlet.

www.amico.com

Amico Pipeline | www.amico.com

85 Fulton Way, Richmond Hill Ontario, L4B 2N4, Canada

71 East Industry Court, Deer Park NY 11729, U.S.A

Toll Free Tel:1.877.462.6426Toll Free Fax:1.866.440.4986Tel:905.764.0800Fax:905.764.0862Email:info@amico.com