

Product Overview

GoLift One Hook Carry Bar



Quick Release One Hook Carry Bar:
CB-GO-LIFT-REG-1HK



Fixed One Hook Carry Bar:
CB-GO-LIFT-REG-1HK-FIX




GoLift One Hook Carry Bar

The one hook carry bar is a mounting bar between the lift and the sling that supports the patient. The carry bar has a maximum capacity of 700 lbs and has one hook on either side to mount the strap loops of a sling. The hooks feature antibacterial silicone latches that prevent the strap loops from accidentally coming off when there is slack on the sling. The steel construction of the carry bar allows for convenient docking of the magnetic hand control when not in use. The entire carry bar is designed with minimal seams for infection control. Please see GoLift manual for a detailed list of acceptable cleaners to wipe down the surfaces.

The **quick connect version of the one hook carry bar** can be removed and reattached from the lift strap with the press of a button. The locking mechanism has a gravity lock that prevents the button from being pressed when there is any weight on the hook (including the carry bar weight). When removing the carry bar, raise the bar so that slack on the strap allows the button to be pressed in.



The **fixed version of the one hook carry bar** can only be used with the GoLift Loop Strap. The strap of the lift tapers into a one inch loop through which a pin can be inserted to secure the carry bar in place. A key ring is inserted through the end of the pin to keep it from sliding out of place.



Variations:

CB-GO-LIFT-REG-1HK:

Quick Release One Hook Carry Bar

Compatible with GoLift systems with a hook at the end of the strap.

CB-GO-LIFT-REG-1HK-FIX:

Fixed One Hook Carry Bar

Compatible with GoLift systems with a loop at the end of the strap

(GO-LIFT-400-LS/ GO-LIFT-700-LS) and cannot be used with the quick connect version.


Features:

Maximum capacity: 700 lbs


Span: 19.175"

Dimensions

Quick Release One Hook Carry Bar



Fixed One Hook Carry Bar



*Dimensions are shown in inches



AMO-PO-GOLIFT-1HOOK-CARRY-BAR 06.22.2022

Amico Mobility Solutions Corporation
122-B East Beaver Creek Road, Richmond Hill, ON L4B 1G6, Canada | 600 Prime Place, Hauppauge, NY 11788, USA
Toll Free Tel: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 | Tel: 905.764.0800 | Fax: 905.764.0862 | www.amico.com

