



# CERTIFICATE

**No. U8 098792 0002 Rev. 00**

**Holder of Certificate:** Amico Accessories Inc.

122 East Beaver Creek  
Richmond Hill ON L4B 4V1  
CANADA

**Certification Mark:**



**Product:** General Medical Devices  
Equipment Mounting Arm

This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL and a Standards Council of Canada accredited Certification body.

**Test report no.:** 7169006435-000

Date, 2020-01-03

( Peter Keith )

Peter Kato



America

# CERTIFICATE

No. U8 098792 0002 Rev. 00

**Model(s):**

AHM-xxxx-wyyz-(a)-vv  
AHM-FITx-bcd-ee-ff  
LIT-www-aa(xx)v-(yy)(zz)-bb  
SSM-xxx-yyz(a)-vv  
KEYBRDw-xxxxyy-(a)  
(See attachment for naming convention)

**Tested  
according to:**

CAN/CSA-C22.2 No. 60601-1:2014 excluding electromagnetic  
compatibility (Clause 17)  
ANSI/AAMI ES60601-1:2005/A1:2012-08 excluding electromagnetic  
compatibility (Clause 17)

**Production  
Facility(ies):  
Parameters:**

098792

Rated Input Voltage: 5Vdc  
Rated Input Current: 1A







America

# CERTIFICATE

No. U8 098792 0002 Rev. 00

AHMM-FITx-bcd-ee-ff – Product Line ("Adjustable Height Monitor – Falcon IT")

Where:

x, Color type, one of:

W, white; or

X, custom, not affecting safety.

b, Keyboard tray, one of:

F, falcon; or

M, mini; or

D, sidekick secondary surface; or

X, custom.

c, Extension length, one of:

0, no extension (D = 5); or

7, 6" (15.24cm) extension; or

3, 9" (22.86cm) extension; or

1, 12" (30.48cm) extension; or

5, 16" (40.64cm) extension.

d, Extension type, one of:

F, folding; or

R, reverse folding; or

5, non-folding extension (C = 1, 3, 5, 7) or no extension (C = 0).

ee, Spring strength to support monitor weight, one of:

25, 5.5-8lbs (2.5-3.6kg); or

35, 7.5-10lbs (3.4-4.5kg); or

50, 13-15.5lbs (5.8-7kg); or

60, 15.5-18lbs (7-8.1kg); or

65, 15.5-18lbs (7-8.1kg); or

75, 18-22lbs (8.1-10kg); or

FH, Fixed height support rod (>10kg, subject to overall loading restriction)

f, Arm mount option, one of:

(For any cd):

A, Amico Horizontal adapter; or

AB, Amico horizontal direct mount - bottom; or

AT, Amico horizontal direct mount - top; or

M or MON or (blank), standard channel adapter; or

D, Dovetail/Ohmeda channel adapter; or

V, Amico vertical rail adapter.

(For cd = 05 only):

S, HSI Mount; or

G, Hill-Rom vertical rail adapter; or

C or COW, CONMED mount.



America

# CERTIFICATE

No. U8 098792 0002 Rev. 00

LIT-www-aa(xx)v-(yy)(zz)-bb - Product Line ("Falcon Lite IT Station")

Where:

**www**, Device that is mounted, one of:

LCD, LCD, VESA brackets provided (aa, yy and v are (blank)); or  
(blank), Monitor shelf or IT station configuration (Either aa and yy are not (blank), or v not (blank)).

**aa**, Accessory tray, one of:

WD, Sidekick dual surface; or  
WF, Falcon keyboard tray; or  
WM, mini-keyboard tray; or  
(blank), LCD, VESA bracket provided option (Either www = LCD or V not (blank)).

**xx**, Arm length and configuration, one of:

45, no extension arm (flush); or  
75, 6" (15.24 cm) single extension arm; or  
25, 9" (22.86 cm) single extension arm; or  
05, 12" (30.5cm) single extension arm; or  
55, 16" (40.64 cm) single extension arm; or  
35, 9" x 9" non-folding dual extension arm; or  
15, 12" x 12" non-folding dual extension arm; or  
65, 16" x 16" non-folding dual extension arm; or  
3F, 9" x 9" dual folding extension arm; or  
1F, 12" x 12" dual folding extension arm; or  
5F, 16" x 16" dual folding extension arm; or  
3R, 9" x 9" Reverse folding extension arm, vertically reversed; or  
1R, 12" x 12" Reverse folding extension arm, vertically reversed; or  
5R, 16" x 16" Reverse folding extension arm, vertically reversed.

**vv**, Shelf type for monitor mounting, (www = (blank), aa = (blank), and yy = (blank))

0, 5" (12.7cm) shelf with 30° tilt limit; or  
2, 6" (15.2cm) drop on shelf; or  
3, Philips IntelliVue Shelf; or  
4, Accessory add on top; or  
6, 5" (12.7cm) shelf for back mounts; or  
7, 5" (12.7cm) shelf with slide-in VESA 75/100 plate; or  
8, 3" shelf with 30° tilt limit; or  
9, 3" shelf; or

W, Welch Allyn Head with 30° tilt limit; or  
(blank), LCD or IT station (Either www not (blank), or aa and yy not (blank)).

**yy**, Spring strength for support of monitor and peripherals, one of:

25, 5.5-8lb (2.5-3.6kg); or  
35, 7.5-10lb (3.4-4.5kg); or  
50, 13-15.5lb (5.8-7.0kg); or  
60, 15.5-18lb (7-8.1kg); or  
65, 15.5-18lb (7-8.1kg); or  
75, 18-22lb (8.1-10kg); or  
FH, Fixed height support rod (>10kg, subject to overall loading restriction); or  
(blank), LCD, or Monitor shelf option (Either WWW = LCD or V not (blank)).

**zz**, Arm mount option, one of:

(For any xx)

A, Amico Horizontal adapter; or  
AB, Amico horizontal direct mount - bottom; or  
AT, Amico horizontal direct mount - top; or  
M or MON or (blank), standard channel adapter; or  
D, Dovetail/Ohmeda channel adapter; or  
V, Amico vertical rail adapter.

(For xx = 45, -5, 25, 55, or 75 only:)

ARS, Amico rail system adapter; or  
S, HSI Mount; or  
G, Hill-Rom vertical rail adapter; or  
C or COW, CONMED mount.

**bb**, Normal or Heavy duty option, one of:

01, Heavy Duty; or  
(blank), Regular.





America

# CERTIFICATE

No. U8 098792 0002 Rev. 00

**SSM-xxx-yyz-(a)-vv** – Product line ("Side to side motion")

Where:

**xxx**, Device that is mounted, one of:

LCD, LCD, VESA brackets provided (z (blank)); or  
(blank), Monitor, shelf provided (z not (blank)).

**yy**, Arm length and configuration, one of:

25, single 9" (22.9 cm) SSM Arm; or  
35, dual 9" (22.9 cm) SSM Arm; or  
05, single 12" (30.5 cm) SSM Arm; or  
15, dual 12" (30.5 cm) SSM Arm; or  
45, no arm (flush mount).

**z**, Shelf type for monitor mounting (xxx (blank)), one of:

0, 5" (12.7cm) shelf with 30° tilt limit; or  
2, 6" (15.2cm) drop on shelf; or  
3, Philips IntelliVue Shelf; or  
4, Accessory add on top; or  
6, 5" (12.7cm) shelf for back mounts; or  
7, 5" (12.7cm) shelf with slide-in VESA 75/100 plate; or  
8, 3" shelf with 30° tilt limit; or  
9, 3" shelf; or  
W, Welch Allyn Head with 30° tilt limit; or  
(blank), not for monitor mounting (xxx not (blank)).

**a**, arm mount option, one of:

(For any yy):

A, Amico Horizontal adapter; or  
M or MON or (blank), standard channel adapter; or  
D, Dovetail/Ohmeda channel adapter; or  
V, Amico vertical rail adapter.

(For yy = 05, 25, 45 only)

F, Fairfield rail adapter; or  
R2, 2" post adapter; or  
E25, 25mm European rail adapter; or  
E35, 35mm European rail adapter; or  
PLV, 0.5"-1.5" Adjustable vertical post adapter; or  
PLH, 0.5"-1.5" Adjustable horizontal post adapter; or  
ARS, Amico rail system adapter; or  
S, HSI Mount; or  
G, Hill-Rom vertical rail adapter; or  
C or COW, Ergotron and Rubbermaid Rail mount.

**vv**, Normal or Heavy duty option, one of:

01, Heavy Duty; or  
(blank), Regular.



America

# CERTIFICATE

No. U8 098792 0002 Rev. 00

## KEYBRDw-xxxxy-(a) - Product Line ("Keyboard Arm")

Where:

w, Keyboard option, one of:

- 4, Falcon with ultra flush; or
- A, Regular Keyboard; or
- D, Dual surface with keyboard; or
- E, Mini keyboard with ultra flush; or
- F, Amico Falcon Keyboard; or
- K, Large Keyboard with ultra flush; or
- L, Large Keyboard; or
- N, Falcon keyboard with flush; or
- M, Amico Mini Keyboard.

xxx, Arm option, one of:

- AHM, Adjustable height arm (AHM); or
- SSM, Side to side monitor arm (SSM); or
- LIT, Falcon Lite IT station (LIT).

yy, Extension option, one of:

(If xxx = AHM only)

- 05, no extension;
- 15, 12" (30.5 cm) single, non-folding extension arm;
- 7F, 6" (15.24 cm) folding extension arm; or
- 3F, 9" (22.86 cm) folding extension arm; or
- 1F, 12" (30.48 cm) folding extension arm; or
- 5F, 16" (40.64 cm) folding extension arm; or
- 7R, 6" (15.24 cm) reverse folding extension arm; or
- 3R, 9" (22.86 cm) reverse folding extension arm; or
- 1R, 12" (30.48 cm) reverse folding extension arm; or
- 5R, 16" (40.64 cm) reverse folding extension arm.

(if xxx = AHM or LIT):

- 75, 6" (15.24 cm) single, non-folding extension arm; or
- 25, 9" (22.86 cm) single, non-folding extension arm; or
- 55, 16" (40.64 cm) single, non-folding extension arm.

a, Arm mount option, one of:

(For xxx = SSM and yy = 2y, 0y, or 4y only):

- S, HSI Mount; or
- G, Hill-Rom vertical rail adapter; or
- C or COW, CONMED mount; or
- F, Fairfield rail adapter; or
- R2, 2" post adapter; or
- E25, 25mm European rail adapter; or
- E35, 35mm European rail adapter; or
- PLV, 0.5"-1.5" Adjustable vertical post adapter; or
- PLH, 0.5"-1.5" Adjustable horizontal post adapter; or
- ARS, Amico rail system adapter.

(For xxx = LIT and yy = 45, 05, 25, 55, or 75 only)

- S, HSI Mount;
- G, Hill-Rom vertical rail adapter;
- C or COW, CONMED mount;
- AT, Amico horizontal direct mount - top; or
- AB, Amico horizontal direct mount - bottom.
- ARS, Amico rail system adapter.

(if xxx = LIT only)

- 45, no extension; or
- 05, 12" (30.5 cm) single, non-folding extension arm; or
- 35, 9" x 9" non-folding dual extension arm; or
- 15, 12" x 12" non-folding dual extension arm; or
- 65, 16" x 16" non-folding dual extension arm; or
- 3F, 9" x 9" dual folding extension arm; or
- 1F, 12" x 12" dual folding extension arm; or
- 5F, 16" x 16" dual folding extension arm; or
- 3R, 9" x 9" dual reverse folding extension arm; or
- 1R, 12" x 12" dual reverse folding extension arm; or
- 5R, 16" x 16" dual reverse folding extension arm.

(if xxx = SSM only)

- 2, single 9" (22.9 cm) SSM Arm; or
- 0, single 12" (30.5 cm) SSM Arm; or
- 3, dual 9" (22.9 cm) SSM Arm; or
- 1, dual 12" (30.5 cm) SSM Arm; or
- 4, no arm (flush mount).

(For any xxx and yy)

- A, Amico Horizontal adapter; or
  - M or MON or (blank), standard channel adapter;
  - D, Dovetail/Ohmeda channel adapter; or
  - V, Amico vertical rail adapter.
- (For xxx = AHM and yy = 05 only)
- S, HSI Mount; or
  - G, Hill-Rom vertical rail adapter; or
  - C or COW, CONMED mount; or
  - AT, Amico horizontal direct mount - top; or
  - AB, Amico horizontal direct mount - bottom.



America

# CERTIFICATE

No. U8 098792 0002 Rev. 00

**MON-DUAL** - Add-on ("Dual Monitor") for installation with two **SSM** bodies, and optimally, one **AHM** body)

Where:

xx, First Extension type, one of:

May optionally be installed with AHM body:

4L, flush mount (no arm)

2L, single 9" SSM arm

Cannot be installed with AHM body:

0L, single 12" SSM arm

yy, Second Extension type, one of:

May optionally be installed with AHM body:

4L, flush mount (no arm)

2L, single 9" SSM arm

Cannot be installed with AHM body:

0L, single 12" SSM arm