

# Capsule™

## Mounting Solutions



Neuron®



Neuron® 2



Axon®



Smartlinx Vital Plus

### Mounting options for every application!



R2  
2" (5.1 cm)  
Post Mount



G  
Hill-Rom Vertical  
Rail Mount



PLV  
0.5"-1.5" (1.3 cm - 3.8 cm)  
Adjustable  
Vertical Post Mount



V  
VRS  
Mount



COW  
Ergotron and Rubbermaid  
Rail Mount



M  
Standard Channel  
Mount



S  
HSI  
Mount



ARS  
Amico Rail System  
Mount



F  
FTA Rail  
Mount



PLH  
0.5"-1.5" (1.3 cm - 3.8 cm)  
Adjustable Horizontal Rail/  
Post Mount



A  
ARS Channel  
Mount



PC1  
0.5"-1.5" (1.3 cm - 3.8 cm)  
Adjustable Rotating Rail/  
Post Mount

Anesthesia Machines · Headwalls · Mobile Computer Carts · Roll Stands · IV Poles



## Carescape V100 Roll Stand Extender

PN: CS-RS-LCD-CLMP-78-12-02

Use to mount on:

GE Carescape V100 or Dinamap Roll Stands as shown in picture

### SPECIFICATIONS:

*Dimensions (without dock or Neuron/Neuron 2)*

Overall: 31.5 cm (12.4") x 11.4 cm (4.5") x 13.5 cm (5.3")

Outer Post Diameter: 3.2 cm (1.25")

Shipping Carton: 33 cm (13") x 18 cm (7") x 18 cm (7")

Net Weight: 0.95 kg (2.3 lbs)

Inner Post Diameter: 2.5 cm (1")

Shipping Weight: 1.4 kg (3 lbs)



## Accutorr Plus Roll Stand Extender

PN: CS-DS-ACC-RS-CLMP

Use to mount on:

Mindray Datascope Accutorr Roll Stands as shown in picture

### SPECIFICATIONS:

*Dimensions (without dock or Neuron/Neuron 2)*

Overall: 21.6 cm (8.5") x 11.4 cm (4.5") x 13.5 cm (5.3")

Outer Post Diameter: 3.2 cm (1.25")

Shipping Carton: 33 cm (13") x 18 cm (7") x 18 cm (7")

Net Weight: 1 kg (2.2 lbs)

Inner Post Diameter: 2.5 cm (1")

Shipping Weight: 1.4 kg (3 lbs)



## Accutorr (Plus or V) Value Stand Extender

PN: CS-DS-ACC-RSV-CLMP

Use to mount on:

Mindray Value Stands as shown in picture

### SPECIFICATIONS:

*Dimensions (without dock or Neuron/Neuron 2)*

Overall: 53.3 cm (21") x 11.4 cm (4.5") x 21.6 cm (8.5")

Outer Post Diameter: N/A

Shipping Carton: 61 cm (24") x 28.6 cm (11.25") x 15.2 cm (6")

Net Weight: 1.9 kg (4.2 lbs)

Inner Post Diameter: N/A

Shipping Weight: 3 kg (6.6 lbs)



## Accutorr (Plus or V) Transport Stand Extender

PN: CS-DS-ACC-RST-CLMP

Use to mount on:

Mindray Transport Stands as shown in picture

### SPECIFICATIONS:

*Dimensions (without dock or Neuron/Neuron 2)*

Overall: 54.1 cm (21") x 12.9 cm (5.1") x 24.9 cm (9.8")

Outer Post Diameter: N/A

Shipping Carton: 71 cm (28") x 25.4 cm (10") x 19 cm (7.5")

Net Weight: 1.9 kg (4 lbs)

Inner Post Diameter: N/A

Shipping Weight: 3 kg (7 lbs)



## VS3 Roll Stand Extender

PN: CS-PH-VS3-RS-CLMP

Use to mount on:

Phillips VS3 Monitor Roll Stands as shown in picture

### SPECIFICATIONS:

*Dimensions (without dock or Neuron/Neuron 2)*

Overall: 52.1 cm (20.5") x 11.4 cm (4.5") x 17.8 cm (7")

Outer Post Diameter: 5.1 cm (2")

Shipping Carton: 61 cm (24") x 28 cm (11") x 18 cm (7")

Net Weight: 2.1 kg (4.6 lbs)

Inner Post Diameter: N/A

Shipping Weight: 3.2 kg (7 lbs)



## LXi Roll Stand Extender

PN: CS-WA-LXI-RS-CLMP

Use to mount on:

Welch Allyn LXi Roll Stands as shown in picture

### SPECIFICATIONS:

*Dimensions (without dock or Neuron/Neuron 2)*

Overall: 53.3 cm (21") x 11.4 cm (4.5") x 17.3 cm (6.8")

Outer Post Diameter: 3.2 cm (1.25")

Shipping Carton: 61 cm (24") x 28 cm (11") x 18 cm (7")

Net Weight: 2.1 kg (4.6 lbs)

Inner Post Diameter: N/A

Shipping Weight: 3.2 kg (7 lbs)



## Connex Roll Stand Extender

PN: CS-WA-SSM2-RS-CLMP

Use to mount on:

Connex Roll Stand Extender as shown in picture

### SPECIFICATIONS:

*Dimensions (without dock or Neuron/Neuron 2)*

Overall: 58.2 cm (22.9") x 11.4 cm (4.5") x 28.7 cm (11.3")

Outer Post Diameter: 3.2 cm (1.25")

Shipping Carton: 71 cm (28") x 25.4 cm (10") x 19 cm (7.5")

Net Weight: 1.9 kg (4 lbs)

Inner Post Diameter: N/A

Shipping Weight: 3 kg (7 lbs)



## Wall Mount

PN: SSM-LCD-15

Use to mount on:

Wall channel or any GCX compatible wall channel

### SPECIFICATIONS

Reach: 61 cm (24")

Net Weight: 2.2 kg (5 lbs)

Shipping Weight: 2.7 kg (6 lbs)

Shipping Carton: 44 cm (17") x 28 cm (11") x 16 cm (6")



## Wall Mount

PN: SSM-LCD-05

Use to mount on:

Wall channel or any GCX compatible wall channel

### SPECIFICATIONS

Reach: 38 cm (15")

Net Weight: 2 kg (5 lbs)

Shipping Weight: 2.3 kg (5 lbs)

Shipping Carton: 44 cm (17") x 28 cm (11") x 16 cm (6")



## Wall Mount

PN: SSM-LCD-45

Use to mount on:

Wall channel or any GCX compatible wall channel

### SPECIFICATIONS

Reach: 11.4 cm (4.5")

Net Weight: 1.1 kg (2.4 lbs)

Shipping Weight: 1.4 kg (3 lbs)

Shipping Carton: 33 cm (13") x 18 cm (7") x 18 cm (7")



## Ultra Flush Wall Mount

PN: PM5-MON

Use to mount on:

Wall channel or any GCX compatible wall channel

### SPECIFICATIONS

Reach: 3.8 cm (1.5")

Net Weight: 0.8 kg (1.8 lbs)

Shipping Weight: 1 kg (2.2 lbs)

Shipping Carton: 22.9 cm (9") x 22.9 cm (9") x 22.9 cm (9")



## Wall Channel

PN: MRS-FL-010

Use to mount on:

For mounting wall mounted capsule solutions;  
SSM-LCD-15, SSM-LCD-05, SSM-LCD-45

### SPECIFICATIONS

Net Weight: 0.3 kg (0.7 lbs)

Shipping Weight: 0.9 kg (2 lbs)

Shipping Carton: 33 cm (13") x 18 cm (7") x 18 cm (7")



## Power Supply Holder

PN: PWRHLD-PC1

Use to mount on:


Poles - 0.5"-1.5" / 12.7 mm - 38 mm / 1.27 cm - 3.8 cm  
DIN Rails - 25 mm / 30 mm / 35 mm

### SPECIFICATIONS

Net Weight: 0.4 kg (1 lbs)

Shipping Weight: 0.9 kg (2 lbs)

Shipping Carton: 16 cm (6") x 16 cm (6") x 11 cm (4")



## Ultra Flush Pole Mount

PN: PM5-PLV

Use to mount on:

Poles - 0.5" - 1.5" / 12.7 mm - 38 mm / 1.27 cm - 3.8 cm  
DIN Rails - 25 mm / 30 mm / 35 mm

### SPECIFICATIONS

Reach: 3.8 cm (1.5")

Shipping Weight: 1 kg (2.2 lbs)

Net Weight: 0.8 kg (1.8 lbs)

Shipping Carton: 22.9 cm (9") x 22.9 cm (9") x 22.9 cm (9")



## DIN Rail and Pole Mount

PN: SSM-LCD-45-PC1

Use to mount on:

Poles - 0.5"-1.5" / 12.7 mm - 38 mm / 1.27 cm - 3.8 cm  
DIN Rails - 25 mm / 30 mm / 35 mm Rails

### SPECIFICATIONS

Net Weight: 1.3 kg (2.9 lbs)

Shipping Weight: 1.8 kg (4 lbs)

Shipping Carton: 33 cm (13") x 18 cm (7") x 18 cm (7")



## Tamper Resistant Mounting Screws

PN: LCD-TPHWR-M4-12

Use to mount on:

Neuron docking station to mount with tamper resistant screws.

Kit includes: (4) M4 tamper resistant screws and specialty screw driver



## Custom Applications

The Capsule Neuron/Neuron 2 is compatible with some of Amico Accessories Inc.'s product lines including the AHM and Roll Stands.



## Ordering Information:

Contact our Global Representatives for more information or for a unique solution.

Amico Accessories Inc. | [www.amico.com](http://www.amico.com)

Tel: +01 905 763 7778

Email: [global@amico.com](mailto:global@amico.com)

