



## Features:

- Load capacity up to 25 lbs (11 kg)
- Overall height 42" (106.7 cm)
- 3" (7.6 cm) locking casters
- 1 or 2 ft (30.4 cm or 60.9 cm) horizontal rail options available
- Custom color options available (white is standard)
- Sturdy 24" (61 cm) base
- Variety of accessories available



Counterweights available for monitor over 5 lbs. Please call to confirm which style fits your need.

### ARS Equipment Roll Stand



- Amico Horizontal Equipment Rail
- Available in 1 or 2 ft (30.4 cm or 60.9 cm) lengths
- 1 ft (30.4 cm) shown

1 Foot - Part Number: ARS-RS042-012  
2 Feet - Part Number: ARS-RS042-024

### FRS Equipment Roll Stand



- Fairfield compatible Equipment Rail
- Available in 1 or 2 ft (30.4 cm or 60.9 cm) lengths
- 1 ft (30.4 cm) shown

1 Foot - Part Number: FRS-RS042-012  
2 Feet - Part Number: FRS-RS042-024

### Equipment Roll Stand



- Using a pole clamp, this roll stand can be completely customized to fit your needs

Part Number: RS5-42

# Roll Stand Accessories



## Hose Assemblies



- Used to connect a medical gas block to a source
- Custom made to order

Call for part numbers

## Gas Blocks



- Snap-on outlets 1, 2 or 3
- Available for all gases
- All outlet fittings available

Amico gas block options shown

## Canisters & Canister Holders



- Disposable and reusable canisters available
- Ring bracket and canister sold separately
- 1200cc, 2000cc, and 3000cc size available

Call for part numbers

## Vacuum Slides



- Conveniently mount your canister suction ring bracket

Part Number: SLIDE-R125

## Blender Bracket



- Used for mounting oxygen/air blenders

Part Number: BLBRKT-R125


## Cord Wrap



- Ideal for wrapping wires, hoses and tubing

Part Number: CORDWRAP-R125

## Wire Basket



- 10" W x 6" D x 9" H white wire basket (25.4 cm x 15.24 cm x 22.86 cm)
- Convenient storage

Part Number: 100609-R125-BB

## Humidifier Bracket



- Used for mounting heated humidifiers

Part Number: HUMBRKT-R125

Please visit [www.amico.com](http://www.amico.com) for all technical specifications, installation instructions, detailed AutoCAD drawings, and more Mounting Solutions.



Distributed By:



Amico Accessories Inc. 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada  
Toll Free Tel: 1.877.264.2697 | Tel: 905.763.7778 | Fax: 905.763.8587 | Email: [accessories@amico.com](mailto:accessories@amico.com)

[www.amico.com](http://www.amico.com)