






























Device	Mounting Kit	Mounting Platform	Restrictions	Add-ons
 <p>GE Solar 8000 Monitor 15", 17", & 19" 12.3 lbs (5.6 kg)*</p>	<p>N/A (The Solar 8000 Flat Panel has a standard VESA mounting pattern)</p>	 <p>AHM-057</p>	<p>N/A</p>	 <p>RAC4</p>  <p>Trim Knob Keypad</p>  <p>Solar 8000 Base Unit</p>
	 <p>GE-8000-SM-FH, MON-VESA-ADA 2.9 lbs (1.3 kg)*</p>	 <p>Solar 8000 Base Unit</p>		<p>Requires an additional mount, see Solar 8000 Base Unit details below</p>

GE Healthcare Ordering Guide Reference Chart


Device	Mounting Kit	Mounting Platform	Restrictions	Add-ons
 <p>Trim Knob Keypad 1 lb (0.5 kg)*</p>	 <p>GE-KEY8K-PLT-00 1.5 lbs (0.7 kg)*</p>	 <p>AHM-057</p>  <p>SSM-057</p>	This unit sandwiches between the flat panel and the AHM/SSM or the GE-8000 SM-FH	N/A
 <p>Solar 8000 Base Unit 12.3 lbs (5.6 kg)*</p>	<p>N/A (The Solar 8000 Base Unit has a mounting plate build into the bottom)</p>	 <p>GE-8000-MON-00</p>		
 <p>RAC4 7.8 lbs (3.4 kg)*</p>		 <p>AHM-050</p>		
		 <p>SSM-050</p>		
	 <p>GE-RAC4-SI-00 1.5 lbs (0.7 kg)*</p>	 <p>AHM-057</p>  <p>SSM-057</p>	Mounts below with primary device on AHM/SSM	N/A

GE Healthcare Ordering Guide Reference Chart

Device	Mounting Kit	Mounting Platform	Restrictions	Add-ons
 <p>Dual RAC4 15.6 lbs (6.8 kg)*</p>	 <p>GE-RAC4-SS-00 2.8 lbs (1.3 kg)*</p>	 <p>AHM-057</p>	Mounts below with primary device on AHM	N/A
 <p>Cuffs and Patient Leads</p>	N/A	 <p>100609-MON-BB</p>	N/A	N/A
 <p>GE Unity 2.4 lbs (1.1 kg)*</p>	 <p>GE-UNITY-PLT-00 0.3 lbs (0.1 kg)*</p>	 <p>GE-8000-MON-00</p>	N/A	N/A
 <p>GE PDM 2.8 lbs (1.3 kg)*</p>	N/A	 <p>GE-PRNPD-LBT-MON</p>	N/A	N/A

GE Healthcare Ordering Guide Reference Chart



Device	Mounting Kit	Mounting Platform	Restrictions	Add-ons
 <p>PRN50 5.7 lbs (2.6 kg)*</p>	 <p>GE-PRNPD-LBT-00 0.4 lbs (0.2 kg)*</p>	 <p>GE-RAC4-SI-00</p>	N/A	N/A
		 <p>GE-RAC4-SS-00</p>		

*Weight provided to help approximate mounting capacity. For reference only.