























GE Healthcare Ordering Guide Reference Chart








Device	Mounting Kit	Mounting Platform	Restrictions	Add-ons
 <p>Transport Pro Monitor 8.4 lbs (3.8 kg)*</p>	 <p>GE-TRPRO-BKT-00 (Allows GE PDM connection to the back of the monitor) 0.5 lbs (0.2 kg)*</p>	 <p>AHM-057</p>	N/A	 <p>RAC2</p>  <p>Power Supply</p>  <p>Cuffs and Patient Leads</p>  <p>GE Unity</p>  <p>GE PDM</p>
		 <p>SSM-057</p>	N/A	
		 <p>TTM-16-457-A00</p>	N/A	
	<p>N/A (The Transport Pro has a 100x100 VESA pattern on the back, but the RAC2 and the GE PDM will not be able to mount to the back of the Transport Pro)</p>	 <p>AHM-057</p>	N/A	
		 <p>SSM-057</p>	N/A	
		 <p>TTM-16-457-A00</p>	N/A	

Device	Mounting Kit	Mounting Platform	Restrictions	Add-ons
 <p>RAC2 5 lbs (2.3 kg)*</p>	N/A	 <p>GE-RAC2M-MON-00</p>	Mounts on MRS rail with AHM/SSM	N/A
	 <p>GE-RAC4-SI-00 1.5 lbs (0.7 kg)*</p>	 <p>AHM-057</p>	Mounts below with primary device on AHM/SSM	
		 <p>SSM-057</p>		
 <p>Power Supply 1.5 lbs (0.9 kg)*</p>	N/A	 <p>PWRHLD-MON</p>	N/A	N/A
 <p>Cuffs and Patient Leads</p>	N/A	 <p>100609-MON-BB</p>	Mounts on MRS rail with AHM/SSM	N/A

GE Healthcare Ordering Guide Reference Chart



Device	Mounting Kit	Mounting Platform	Restrictions	Add-ons
 <p>GE Unity 2.4 lbs (1.1 kg)*</p>	 <p>GE-UNITY-PLT-00 0.3 lbs (0.1 kg)*</p>	 <p>GE-8000-MON-00</p>	N/A	N/A
 <p>GE PDM 2.8 lbs (1.3 kg)*</p>	<p>N/A (The GE PDM mounts to the back of the Transport Pro only with GE-TRPRO-PLT-00)</p>	N/A	N/A	N/A
	N/A	 <p>GE-PRNPD-LBT-MON</p>	Should not be used with GE-TRPRO-BLT-00	N/A

*Weight provided to help approximate mounting capacity. For reference only.