



## Gooseneck Tablet Roll Stand




Amico

# Telehealth Workstations

## Gooseneck Roll Stand



Developed specifically for healthcare in collaboration with healthcare facilities, Amico's Gooseneck Roll Stand combines utility, durability, and customization into a robust piece of equipment. Easy to use, secure, and easy to clean, Amico's Gooseneck Roll Stand is a Telehealth Communication Station applicable for use by patients and caregivers.




# Standard Equipment



## Locking Tablet Holder

Secure, key-lock tablet holder for tablets 11" - 13.5" (27.94 cm - 34.29 cm) (outside, corner to corner dimension), up to 0.4" (1.01cm) thick  
(See Part Matrix on pg. 5 for factory-installed tablet size options, and Large Tablet Kit on pg. 5 for on-site upgrades)



## Tablet Adjustment

360° tablet adjustment for optimal viewing position, with hand-tight securement



## Gooseneck

Infection-control friendly, wipeable, medical-grade gooseneck



## Reach

23" (58.4 cm) long gooseneck for comfortable viewing by patients at center of bed



## Infection Control

Medical-grade wipeable surfaces on gooseneck, pole, and base. Safe for cleaning with commonly used medical cleaners per infection control standards



## Anti-Tip Base

24.5" (62.2 cm) low center of gravity base, compliant to IEC606061 healthcare stability standard



## Effortless Rolling

Large casters designed for hospital use with ergonomic, easy-to-use, brake levers

# Optional Equipment




## Integrated Power Adapter

Secure power coupling with optional pre-run USB-C or Lightning cable (see pg. 5) and 3' (0.9 m) coiled power cord (extends to 6' (1.8 m))




## Plug Type

3' (0.9 m) power cord with plug type available for all international standards (see part matrix on pg. 5 for choices)



## Small Basket Package

5" x 5" x 4"  
(12.7 cm x 12.7 cm x 10.2 cm)  
basket with cord wrap and handle



## Large Basket Package

10" x 6" x 9"  
(25.4 cm x 15.2 cm x 22.9 cm)  
basket with cord wrap and handle



## Large Tablet Holder


Secure locking holder for 13.5" - 15.5" (34.29 cm - 39.37 cm) (outside, corner-to-corner dimension), up to 0.8" thick (See Part Matrix on pg. 5)

# Accessories




## Base Counterweight

Increases stability when adding accessories and storing supplies on your roll stand




## Ergonomic Handle

Ergonomic, infection control friendly handle makes moving the roll stand easy




## Small Basket

5" x 5" x 4"  
(12.7 cm x 12.7 cm x 10.2 cm)  
sterilizable wire basket. Includes cord wrap




## Large Basket

10" x 6" x 9"  
(25.4 cm x 15.2 cm x 22.9 cm)  
sterilizable wire basket. Includes cord wrap



## Cord Wrap

Neatly store your cables and cords when not in use. Attaches to basket or cart



## Chart Holder

Locate your charts close at hand, behind the worktop



## Glove Box Holder

Storage for boxes of gloves (works with tissue boxes too)



## Sanitizer Holder

Securely store your sanitizer pump



## Equipment Wipes Holder

Mounting bracket for your equipment wipes dispenser




## Power Bar

Provide power to multiple devices on the roll stand - all international power types available (power cord included)



## Power Supply Holder


Strap in storage for your power supply or bar




## Large Tablet Kit

Increases tablet holder size range to 13.5" to 15.5" (34.29 cm to 39.37 cm), up to 0.8" (2.04 cm) thick

# Configure your Goose



# Technical Specifications

<b>Dimensions</b>	24.5" x 24.5" x 65.5" (62.2 cm x 62.2 cm x 164.4 cm)
<b>Weight</b>	21.5 lbs (9.8 kg)
<b>Max Tablet Weight</b>	3 lbs (1.36 kg)
<b>Fits Tablet Sizes</b>	<p>Standard: 11" - 13.5" (27.94 cm - 34.29 cm) (outside, corner-to-corner dimension), up to 0.4" thick</p> <p>Large Tablet Kit: 13.5" - 15.5" (34.29 cm - 39.37 cm) (outside, corner-to-corner dimension), up to 0.8" thick</p> <div>  <p>The diagram consists of two side-by-side rectangles representing the tablet holder. The left rectangle is labeled 'WITH STANDARD TABLET MOUNT' and contains a diagonal line from the bottom-left corner to the top-right corner, with the text '11" - 13" TABLET SIZE' below it. The right rectangle is labeled 'WITH EXTENDED KIT' and contains a similar diagonal line, with the text '13" - 15.5" TABLET SIZE' below it. Both rectangles have a small circle with a crosshair in the top-left corner.</p> </div>
<b>Locking Tablet Holder</b>	Yes
<b>Tablet Adjustment</b>	Tilting, 360-degree rotating
<b>Tablet Adjustment Lock</b>	Hand-tight
<b>Tablet Min / Max Height</b>	25.75" / 65.5" (65.4 cm / 166.4 cm)
<b>Gooseneck Length (Stem to Tablet)</b>	22.5" (60.3 cm)
<b>Gooseneck Reach (Horizontal)</b>	16" (40.6 cm)
<b>Gooseneck Material</b>	Wipeable sealed medical-grade sheath covering ridged steel wire.
<b>Casters</b>	Large locking casters designed for hospital use with ergonomic, easy-to-use, brake levers.
<b>Cable Management</b>	Through Gooseneck to optional integrated tablet power adapter with 3' (0.9 m) compact, coiled medical-grade power cord.
<b>Infection Control</b>	Medical-grade wipeable surfaces on Gooseneck, pole, and base. Safe for cleaning with commonly used medical cleaners per infection control standards.
<b>Standards Compliance</b>	<ul style="list-style-type: none"> <li>• ADA</li> <li>• IEC60601 Stability</li> </ul>
<b>Quality Management</b>	<ul style="list-style-type: none"> <li>• ISO 9001:2015</li> <li>• ISO 13485:2016</li> <li>• CE</li> </ul>
<b>Warranty</b>	1 year

Amico Accessories Inc. | 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada  
Toll Free Tel: 1.877.264.2697 | Tel: 905.763.7778 | Fax: 905.763.8587  
Email: accessories@amico.com | www.amico.com

Represented By: