

Mobile Computer Workstation Hummingbird CTM Laptop Cart



Technical Specifications

Load Capacity

Front Load: 25.0 lbs (11.34 kg)

Rear Load: 5.0 lbs (2.27 kg) with max. 2.0" (5.08 cm) offset from the channel surface

Workflow

Modular and flexible to accommodate changing workflows

Maintenance

- Designed with infection control and prevention standards using non-porous, easy-to-clean materials (aluminum, steel, ABS plastic, and nylon)
- Use Caviwipes or other disinfectant wipes (no abrasives, solvents, or trichloroethylene)

Ergonomics

Optional directional caster to enhance steering ability

Cable Management

- Concealed with an integrated easy-to-clean channel cover and cable management throughout the length of the channel/track
- Covers can be easily removed and reinstalled

Base

Sturdy H-shape low profile metal base

Caster

2 total locking casters and 2 swivel casters

Movement

- Laptop Holder Swivel Range: +/- 60°
- Mouse Tray Left-to-Right Range: +/- 5.9" (15.04 cm)

Key Dimensions

Outside Dimensions: 20.5" (52.02 cm) W x 27.0" (68.65 cm) L x 45.3" (115.02 cm) H

CTM Channel Usable Length: 36.6" (93.07 cm)

Laptop Holder Surface: 14.0" (35.70 cm) W x 18.0" (45.70 cm) L

Mouse Tray Surface (left/right): 5.9" (15.04 cm) W x 8.9" (22.48 cm) L

Base: 20.0" (50.79 cm) W x 20.2" (51.36 cm) L

Footprint: 24.9" (63.15 cm) W x 25.9" (65.66 cm) L

Ground Clearance: 6.3" (16.07 cm)

Wheels: 4x 5.0" (12.50 cm) diameter casters

Physical Weight: 60.0 lbs (27.22 kg)

Shipping Dimensions

Box 1 (Column): 5.0" (11.60 cm) W x 5.0" (11.60 cm) L x 36.0" (91.44 cm) H

Box 2 (Base): 22.5" (57.15 cm) W x 22.5" (57.15 cm) L x 10.0" (25.40 cm) H

Box 3 (Laptop Holder): 20.0" (50.80 cm) W x 24.0" (60.96 cm) L x 8.0" (20.32 cm) H

Shipping Weight

Box 1 (Column): 15.0 lbs (6.80 kg)


Box 2 (Base): 35.0 lbs (15.88 kg)

Box 3 (Laptop Holder): 12.0 lbs (5.44 kg)

Warranty


5-Year Warranty

Figure 1: Mouse Tray Adjustments



Drawings

Figure 2: Hummingbird CTM Laptop Cart (48" H) With Flush Side-to-side Laptop Holder (Part No.: HMC-L23W-LP00-00)



Please visit www.amico.com for all technical specifications, installation instructions, detailed AutoCAD drawings, and more Mounting Solutions.



Distributed By:

