



















Philips Ordering Guide Reference Chart



Device	Mounting Kit	Mounting Platform	Restrictions	Add-ons
 <p>Philips IntelliVue MP60/70 22.1 lbs (10 kg)*</p>	 <p>N/A (Philips supplied Table Top Mounts)</p>	 <p>AHM-053</p>	N/A	 <p>FMS Module</p>  <p>Multi-Measurement Module (MMS)</p>  <p>Cable and Patient Leads</p>
		 <p>SSM-053</p>	N/A	
		 <p>TTM-09-453-A00</p>	N/A	

Philips Ordering Guide Reference Chart



Device	Mounting Kit	Mounting Platform	Restrictions	Add-ons
 <p>FMS Module 7.7 lbs (3.5 kg)*</p>	 <p>N/A (Clamp from Philips)</p>	 <p>AHM-053</p>	Clamps below with primary device on AHM/SSM	N/A
		 <p>TTM-09-453-A00</p>	Clamps below with primary device on TTM	N/A
 <p>Multi-Measurement Module (MMS) 1.4 lbs (0.7 kg)*</p>	<p>N/A (Mounts to the back of the MP60/70)</p>	N/A	N/A	N/A
 <p>Cable and Patient Leads</p>	 <p>CRDHLDR2</p>	 <p>AHM-053</p>	Mounts below with primary device on AHM/SSM	N/A
		 <p>SSM-053</p>		
	 <p>101004-MON-BB</p>	N/A	Mounts on MRS rail with AHM/SSM	N/A

*Weight provided to help approximate mounting capacity. For reference only.