















# Welch Allyn Ordering Guide Reference Chart












Device	Mounting Kit	Mounting Platform	Restrictions	Add-ons
 <p><b>Welch Allyn Connex Spot</b>            Standard: 2.9 lbs (1.3 kg)*            Baun: 3.7 lbs (1.7 kg)*            SureTemp: 3.5 lbs (1.6 kg)*</p>	<p>N/A            (The Connex Spot has a standard VESA mounting pattern)</p>	 <p>AHM-057</p>	<p>N/A</p>	 <p>Cable and Patient Leads</p>  <p>7000-PS Power Supply</p>
		 <p>SSM-057</p>	<p>N/A</p>	
		 <p>TTM-08-457-A00</p>	<p>N/A</p>	
		 <p>RS5-427-PKG</p>	<p>N/A</p>	

# Welch Allyn Ordering Guide Reference Chart



Device	Mounting Kit	Mounting Platform	Restrictions	Add-ons
 <p><b>7000-PS (Power Supply)</b></p>	 <p>PWRHLD-MON</p>	 <p>AHM-057</p>	<p>Mounts on MRS rail with AHM/SSM</p>	<p>N/A</p>
		 <p>SSM-057</p>		
	 <p>PWRHLD-R125</p>	 <p>RS5-427-PKG</p>	<p>Mounts on the post of Roll Stand/Table Top Mount</p>	<p>N/A</p>
		 <p>TTM-08-457-A00</p>		

Device	Mounting Kit	Mounting Platform	Restrictions	Add-ons
 <p><b>Cable and Patient Leads</b></p>	 <p>CRDHLDR</p>	 <p>AHM-057</p>	Mounts below with primary device on AHM/SSM	N/A
		 <p>SSM-057</p>		
	 <p>101004-MON-BB**</p>	 <p>AHM-057</p>	Mounts on MRS rail with AHM/SSM***	N/A
		 <p>SSM-057</p>		
	 <p>100609-R125-BB</p>	 <p>RS5-427-PKG</p>	Mounts on Roll Stand****	N/A

\*Weight provided to help approximate mounting capacity. For reference only.

\*\*Or other applicable sizes storage baskets.

\*\*\* Can also mount on RST using C clamp.

\*\*\*\* One Basket with CordWrap on the side is included in RS5-427-PKG.