## ICU Series B-E-AC802 Bed





### **Key Features:**

- 550 lbs (250 kg) weight capacity
- Standard 5th wheel for improved mobility
- 21" (53 cm) low height
- Storage area in base for two oxygen E-cylinder tanks and belongings
- Radiolucent fowler section
- Dual brake pedal mechanisms
- Adjustable knee gatch for enhanced comfort
- Fold-down siderails
- 8" (20 cm) wheel casters
- Pressure relief surfaces

#### **Standard Features:**

- LINAK<sup>™</sup> DC Actuators
- Modular battery backup
- LED night light
- Four IV pole sockets
- Four integrated roller bumpers
- Metal litter deck
- X-ray cassette holder

- CPR release system
- Trendelenburg and reverse Trendelenburg
- Patient controls on inner siderails
- Flat litter deck for compatibility with standard patient surfaces
- · Adjustable bed height
- Adjustable foot extension

- Adjustable litter with chair position
- Integrated digital scale
- Patient exit alarm
- Scale and fowler angle indicator
   on siderail
- Foley bag holder
- Integrated push handles
- 5th wheel for improved mobility

The new Amico ICU bed features a versatile design suitable for a wide range of patients. The ICU bed is differentiated by its enhanced chair position as well as its ergonomic and robust design. This product truly helps provide a more meaningful level of care to patients.

This bed has a battery back-up system to enable motion even when the bed is disconnected or during patient transport. The ICU bed also features LINAK<sup>™</sup> DC actuators that offer a reliable yet silent motion for the bed's electromechanical system.

# **Technical Specifications**



Nodel Number			
Model Number	B-E-AC802		
Overall Dimensions			
Overall Length			
Retracted	83" (211 cm)		
Extended	91" (231 cm)		
Overall Width			
Siderails Up	42" (107 cm)		
Siderails Down	38" (96.5 cm)		
Patient Surface			
Retracted	32" x 71" (81 cm x 180 cm)		
Extended	32" x 79" (81 cm x 201 cm)		
Height Range (to the surface of the litter)			
High	38" (96.5 cm)		
Low	20.5" (52 cm)		
itter Positioning			
Trendelenburg/Reverse Trendelenburg	± 18°		
Backrest	0° - 70°		
Knee Gatch	0° - 21°		
Veight Capacity & Caster Diameter			
Weight Capacity	550 lbs (250 kg)		
Caster Diameter	6" (15 cm)		
lectronics			
Volts			
Domestic	120 VAC		
International	220 VAC		
Ampere Rating			
Domestic	7.5 A		
International	4.0 A		
Current Leakage	<300 µA		
-			
Frequency	50 - 60 Hz		

Optional Features & Accessories:		
Pressure relief surfaces	Overbed table	Bed exit alarm
<ul> <li>Removable IV pole - 2 stage</li> </ul>	<ul> <li>Brake not set alarm</li> </ul>	<ul> <li>Fowler angle indicator</li> </ul>
<ul> <li>IV hook</li> </ul>	<ul> <li>Communication to nurse call systems</li> </ul>	<ul> <li>Hand pendant controls</li> </ul>
Siderail covers	Oxygen tank holder	• 3-Zone bed exit alarm

Amico Beds Corporation reserves the right to change product specifications without notice.

#### Warranty: One year parts only.

Amico Beds Corporation | 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada | 600 Prime Place, Hauppauge, NY 11788, USA Toll Free Phone: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 | Tel: 905.764.0800 | Fax: 905.764.0862 | www.amico.com AB-SP-B-E-AC802-BED 12.08.2022

