



Key Features & Benefits:

- Completely sewed with no zipper
- Economy firm mattress
- Stretch cover is fluid resistant, fire retardant, antifungal and passes antimicrobial test AATCC 30-2004
- Lightweight mattress design allows for ease of rotation and head to foot reversibility

Technical Specifications

Part Numbers


80" x 35" x 6"

M-THEA-WL001-B80356

84" x 35" x 6"

M-THEA-WL001-B84356

Dimensions*



*Custom sizes are available upon request

Amico Beds Corporation reserves the right to change product specifications without notice.