

Care Series

Manual Patient Transfer Stretcher

S-M-100



Key Features:

- Centrally controlled brake system
- Gas spring supported back rest
- Adjustable back rest angle
- Retractable fifth wheel steering system
- Integrated oxygen bottle holder
- 550 lbs (250 kg) weight capacity

Applications:

- Emergency Room
- Recovery
- Pre-Op
- Post-Op

The Amico Emergency Care Stretcher Collection provides caregivers the control they need in the emergency environment. The stretcher frame is made of high quality aluminum and the cart surface and side rails are aesthetically shaped in plastic.

The height adjustment of the stretcher is controlled through the use of a manual lever located at the foot of the stretcher. The back rest angle can be adjusted by depressing a handle to activate the gas-powered spring.

Technical Specifications



Model Number

Model Number	S-M-100
--------------	---------

Overall Dimensions

Overall Length	83" (211 cm)
Overall Width	28" (71 cm)
Patient Surface	28" x 75" (71 cm - 211 cm)
Height Range (to the surface of the litter)	
High	34" (86 cm)
Low	22" (56 cm)

Litter Positioning

Trendelenburg/Reverse Trendelenburg	n/a
Maximum Fowler Section Inclination	75° + 5°
Maximum Knee Section Inclination	n/a

Weight Capacity & Caster Diameter

Weight Capacity	550 lbs (250 kg)
Caster Diameter	8" (20 cm)

Optional Features & Accessories:

- Overbed tables
- Removable IV pole

Amico Beds Corporation reserves the right to change product specifications without notice.

Warranty: One year parts only.