



## Key Features:

- 550 lbs (250 kg) weight capacity
- 24" (61 cm) low height
- Storage area for patient belongings
- Optional radiolucent fowler section
- Easily accessible brake pedals
- Fold-down siderails
- 8" (20 cm) wheel casters
- Pressure relief surfaces
- Articulating headpiece

The Amico Titan EYE/ENT/ORAL Care Surgery Stretcher is designed with functionality and convenience in mind, allowing caregivers to provide a more meaningful level of care to patients.

This product offers key features such as an articulating headpiece that is a stainless steel design for enhanced motion, durability and sterilization; a deep well storage area for patient belongings and an overall ergonomic design to help improve efficiency of healthcare. Arm support is adjustable on either side of the patient.

The Titan EYE/ENT/ORAL Care Surgery Stretcher can also serve as an interim transport stretcher for moving patients between departments.

## Standard Features:

- Four IV pole holders
- Roller bumpers
- Adjustable fowler section
- Fifth wheel steering
- Fold-down siderails
- 8" (20 cm) wheel casters
- Two side mounted pedals for patient surface height adjustment
- Hydraulic lift system capable of supporting up to 550 lbs (250 kg)
- Central braking system
- Ergonomic push handles
- Side bumpers
- Head and foot-end mattress retainer
- Articulating headpiece
- Wrist rest support
- Adjustable arm support

# Technical Specifications



## Model Number & Approval

Model Number	S-AM1-3004
Agency Approval	cULus

## Overall Dimensions

Overall Length	86" (218 cm)
Overall Width	33" (84 cm)

## Litter Positioning

Trendelenburg/Reverse Trendelenburg	+/-16°
Maximum Fowler Section Inclination	80°
Maximum Knee Section Inclination (optional)	35°

## Height Range

Lowest	24" (61 cm)
Highest	33" (83.8 cm)

## Headpiece

Height Adjustment	+4" (10 cm)
Side to Side Adjustment	±15°

## Surface Dimensions

Surface Width	26" (66 cm)
Surface Length	83" (210 cm)
Surface Thickness	4" (10 cm)

## Weight Capacity & Caster Dimension

Maximum Weight Limit	550 lbs (250 kg)
Caster Diameter	8" (20 cm)

## Optional Features & Accessories:

- 3" (8 cm) ultra comfort mattress
- 4"/5" (10 cm/12.7 cm) high density foam mattress with or without fire barrier
- Foley bag holder
- Disposable waste basins
- Utility tray
- Permanent IV pole

Amico Beds Corporation reserves the right to change product specifications without notice.

**Warranty:** One year parts only.

