

Stretchers Titan Series

The Amico Titan Series Transport Stretcher uses an innovative design and incorporates distinctive features aimed to help caregivers provide a more meaningful level of care to patients.

This product is specifically designed to meet the needs of customers by offering...

- 21" (53 cm) low height (best-in-class)
- Standard fifth wheel for improved mobility
- 8" (20 cm) integrated wheel casters
- Integrated foot push handle
- Highest capacity storage tray in the industry
- Storage area for 2 oxygen E-cylinder tanks in the base cover
- Radiolucent fowler board



With the increased awareness regarding fall prevention, Amico is proud to offer a Transport Stretcher that is the best-in-class low height stretcher in the market.

Features



Titan Series



Head Section Articulation




Knee Surface Articulation Crank



Push Handles



Hi-Lo Function Foot Pedal




5th Wheel



Adjustable Side Rails



Brake/Steer Foot Pedal



Permanent Telescopic IV Pole