

SELECT SERIES SECURITY WALL CONSOLE

DRAWING # 0X


(M/N: W-X-SEC-PCON-03)

TYPE: _____
 QUANTITY: X



IF MIRRORED UNITS ARE REQUIRED THE CUSTOMER MUST INDICATE THE QUANTITY OF UNITS BELOW AND AMICO WILL PROVIDE A SEPARATE DRAWING SHOWING THE MIRRORED LAY OUT AND QUANTITY

QUANTITY MIRRORED: _____

SYSTEM DETAILS			
SYMBOL	QTY.	DESCRIPTION	
SEC	1	STAINLESS STEEL SECURITY DOOR	
SLM	1	VACUUM SLIDE - DOOR MOUNTED	
O	1	AMICO GAS OUTLET - OXYGEN	
V	1	AMICO GAS OUTLET - VACUUM	



- NOTE:
1. PUNCH OUT CENTRE TO CENTRE IS 2 1/2"
 2. THE BACK BOX WILL INCORPORATE 1/2" AND 3/8" KNOCKOUTS ON THE TOP AND BOTTOM. THE BACK BOX IS PROVIDED WITH ADJUSTABLE FLANGES TO ACCOMMODATE VARIOUS WALL THICKNESS.
 3. FASCIA WILL BE EXTRUDED ANODIZED ALUMINUM.



IMPORTANT: PLEASE VERIFY THAT THE ABOVE INFORMATION IS CORRECT, AND PROVIDE THE REMAINING DETAILS.

APPROVAL - PRINT AND SIGNATURE _____

DATE _____

PHONE NO. _____

Amico
 85 Fulton Way,
 Richmond Hill, Ontario
 L4B 2N4, CANADA
 Toll-Free: 1-877-462-6428 (T)
 Tel: (905) 764-0960
 Fax: (905) 764-0962
 www.amico.com

HOSPITAL	HOSPITAL NAME	A. NURSE CALL MFGR.: _____	MODEL NO.: _____	DRWG. NO. 0X
LOCATION	LOCATION	B. MEDICAL GAS MFGR.: AMICO	TYPE CONNECTION: _____	DRAWN BY: XX
QTY. X	UNITS AS SHOWN	C. FINISH: SEE NOTE #3	CEILING HEIGHT: _____	CHECKED BY: _____
				REV. NO. 00
				DATE: MM/DD