



CONSOLE OUTLET QUICK DISCONNECT PURITAN - BENNETT COMPATIBLE




TOP VIEW

Wall thickness
may vary
up to 3/4"
[19mm]




SIDE VIEW


1/2" O.D. (12.7mm) (3/8" nominal) size
type "K" copper pipe, rotates 360 deg.
for entry from any angle



FRONT VIEW
BACK ASSEMBLY



FRONT VIEW
LATCH-VALVE ASSEMBLY



SIDE VIEW
LATCH-VALVE ASSEMBLY

MODEL NUMBERS:

The XXX defines the Gas:

- Oxygen = OXY
- Medical Air = AIR
- MedVac = VAC
- Nitrous Oxide = N2O
- Waste Anesthetic Gas Disposal = WAG (NFPA)

O-PBCON-L-XXX

The L defines the Language:

- U = English (NFPA) Oxy, Air, Vac and WAG only.
- S = Spanish (NFPA) Oxy, Air, Vac and WAG only.

Example:

English NFPA - Puritan-Bennett compatible, Oxygen = O-PBCON-U-OXY

85 Fulton Way
Richmond Hill, Ontario
L4B 2N4, CANADA
Tel: (905) 764-0800
Fax: (905) 764-0862

www.amico.com