


CEILING OUTLET DISS SHOWN

MODEL NUMBERS:

The XXX defines the Gas:

Oxygen	= OXY
Medical Air	= AIR
Nitrous Oxide	= N2O
Nitrogen*	= NIT
Carbon Dioxide**	= CO2
Carbogen*	= CAR
MedVac	= VAC
Waste Anesthetic Gas Disposal	= WAG (NFPA)
* = Available only in DISS	
** = Available only in DISS and Ohmeda	





- DIS = DISS
- QD = OHMEDA compatible
- CH = CHEMETRON compatible
- PB = PURITAN - BENNETT compatible

O-DISCEI-L-XXX

The L defines the Language:
 U = English (NFPA) Oxy, Air, Vac and WAG only
 S = Spanish (NFPA) Oxy, Air, Vac and WAG only

Example:
 English NFPA - DISS, Oxygen = O-DISCEI-U-OXY

1/2" O.D. (12.7mm) (3/8" nominal) size
 type "K" copper pipe



SIDE VIEW

Note:
 Amico recommends that you have a blank space between outlets and slides.
 Maximum 3 outlets may be ganged together without additional support.



85 Fulton Way
 Richmond Hill, Ontario
 L4B 2N4, CANADA
 Tel: (905) 764-0800
 Fax: (905) 764-0862
www.amico.com