




Dual Relief Valves CSA Z7396.1

Flow Chart



NOTE: Dimensions are approximate.

- 1 - 1/2" NPT Threads
- 2 - 1/2" Isolation Valves
- 3 - L 700 Regulator
- 4 - Bleed Valve
- 5 - Coupler
- 6 - Pressure Gauge
- 7 - Line Pressure Gauge
- 8 - Relief valve



Note: Nitrogen requires a 300 psi pressure gauge and all other gases require 100 psi. The 1/2" dual line could be adapted to the 1", 2" or 3" supply line by the way adaptors at the in-let and out-let. (adaptors supplied by others)

Model Numbers

M-DLRS-CSA-05 - 1/2" Dual Line Assembly

M-DLRS-CSA-05HP - 1/2" Dual Line Assembly for High Pressure