

Hose Assembly Configurator

Input the hose length by number of feet (12" = 1ft)	<input type="text"/>
Select the medical gas code requirement	<input type="text"/>
Select the medical gas	<input type="text"/>
Select the inlet connection*	<input type="text"/>
Select the outlet connection*	<input type="text"/>
Select if the hose is conductive	<input type="text"/>
Select the size of the hose**	<input type="text"/>

***Note:** Legend of parts for hose colors as well as inlet and/or outlet connections is included in the following pages of this document. Nitrogen is only available for DISS and Schrader connections. Instrument Air is only available for DISS connections. Carbon Dioxide is NOT available for Puritan Bennett connections.

****Note:** 5/16" hose size is only available for VAC, WAGD and AGSS gases at this time.



Example shown:
S-HS-06UO-DFDFC4

Your Part Number

S-HS- - - - - - -

Signature: _____

Signature above confirms that information selected is correct.

Legend of Parts

Hoses



Oxygen (NFPA)



Oxygen (CSA)



Medical Vac (NFPA)



Medical Vac (CSA)



Medical Air (NFPA)



Medical Air (CSA)



Nitrous Oxide



Nitrous Oxide



Waste Anesthetic Gas Disposal (NFPA)



Anesthesia Gas Scavenging System (CSA)



Nitrogen (NFPA)



Nitrogen (CSA)



Carbon Dioxide (NFPA)



Carbon Dioxide (CSA)



Instrument Air (NFPA)



Instrument Air (CSA)

Legend of Parts

Inlet and Outlet Connections



British Male



1/8" MNPT



Chemetron Female



1/4" MNPT



Chemetron Male



Medstar Male



DISS Male (With Demand Check)



NIST Female



DISS Hex Nut



Ohmeda Female



DISS Handtight



Ohmeda Male



DISS Male (Without Demand Check)



Puritan Female



1/8" FNPT



Puritan Male



1/4" FNPT



Schrader Female



Schrader Male