

Spring Arm Camera

Specifications



HD Camera Specifications:

- Up to 120X stable zoom
- Full HD 1080/6p (1920 x 1080) picture quality
- Analog: Component (Y/PB/PR)
- Weight: 210 g (7.4 oz)
- 3,000,000 pixel resolution CMOS chip
- Working distance from 10 mm to 1000 mm
- 2D/3D Noise Reduction
- Dimensions: 45.6 mm x 48.8 mm x 78.0 mm (1³/₁₆" x 1⁵/₁₆" x 3¹/₈")
- Power Consumption: 20 W at 12VDC
- Electricity (VAC): Universal VAC
- Certificates: TUV, FDA

Features:

- Houses HD camera (EV7100)
- Mounts directly to Amico Spring Arm
- Standalone camera unit, no lighthead required
- Allows 360 degrees rotation of camera
- Eco-friendly and long lasting
- Suspension system

