

# V-Series Ceiling Pendant

## Specifications



### Features:

- Hinged back door provides easy access to all internal components allowing for hassle-free maintenance
- Console is 100% customizable; selections range from types of electrical outlets, medgas services, as well as the number of data outlets
- Electric braking system – ergonomic control handle on the console allows for easy maneuverability and no drift when brakes are engaged
- Integrated Equipment Rails allow for positioning of multiple accessories directly onto eight available equipment rails

### Specifications:

- Vertical Range: +/- 22.5 degrees from horizontal
- Rotation Capability: 340 degree
- Max Pay Load: 330 lbs\* (150 kg)
- Input Voltage: 120 – 240 VAC

*\*Based on the configuration, Payload does not account for service head*

Integrated Equipment Rails



Seamless Infection Control



Adjustable Shelves

