

# Aneroid

Wall, Rail and Mobile

6" (15.2 cm) Luminescent Dial



Large Independent Cuff Basket

Latex-free Inflation System



Part Number: AM-DB-LF2118

## Features

- Luminescent dial is readable in low light
- Accurate to within 3 mmHg
- Large 6" (15.2 cm) dial
- Recessed dial face reduces parallax errors
- Lens minimizes reflection to diminish distracting mirror effect
- Thermally treated diaphragm maintains calibration
- 110° gauge swivel for easy viewing
- Available with latex-free blood pressure cuff and 8 ft (2.4 m) coiled tubing inflation system
- Reliable jeweled movement
- No pin stop movement
- Mercury free
- Lifetime calibration warranty
- Independent basket allows convenient placement for increased usability
- Large basket easily holds multiple cuff sizes
- Countless wall and rail mounting configurations available

Amico Diagnostic Incorporated uses precision manufactured components to construct BP Aneroid Sphygmomanometers, providing an accurate reading within 3 mmHg. The clarity of the dial makes it easy to read, and features a scratch and shatter resistant lens. Available with a wide variety of blood pressure cuff and connector configurations.



Part Number: AM-AL-LF2118

# Mobile Mounting



## Features

- Luminescent dial is readable in low light
- Accurate to within 3mmHg
- Large 6" (15.2 cm) dial
- Recessed dial face reduces parallax errors
- Lens minimizes reflection to diminish distracting mirror effect
- Thermally treated diaphragm maintains calibration
- 110° gauge swivel for easy viewing
- Mercury free and latex-free
- Adjustable height from 34" to 48" (86.4 cm to 121.9 cm)
- 18" (45.7 cm) diameter weighted base for added stability
- Spring loaded post helps prevent equipment damage
- Integrated cuff storage
- Aneroid face angled upwards at 20° for ease of viewing

## Rail Mounting

Add flexibility while simplifying installation by rail mounting to your existing headwalls, flatwalls, equipment booms and rail systems.

### L-Bracket Mount



- Aneroid and basket mounted together on a single bracket
- Reduces congestion on the equipment rails
- Available with a variety of rail adapters

Part Number: AM-AL-LF2118

### Horizontal and Vertical Rail Mounts



Part Number: AM-VB-LF2118



Part Number: AM-AB-LF2118

- Easy to move accessory adapters
- Independent cuff basket allows proper placement for ease of access
- Available with a variety of rail adapters

# Aneroid Anatomy



# Aneroid Accessories

## Cuff Basket



- Low profile cuff storage
- Multiple sizes and mounts available

Call for Part Number

## Stethoscopes



- Stainless steel double-head
- Non-chill rim
- Adult and pediatric available
- Latex-free

Call for Part Number

## Two Piece Cuff



- Color coded nylon cover
- Latex-free bladder
- Variety of sizes available
- Connector option to fit any NIBP device

Call for Part Number

## Cuff Connectors



- Large variety of connectors and tubing

Call for Part Number

## Coil Tubing



- Latex-free coil tubing
- Expands to 8 ft (2.4 m)
- Optional connectors

Part Number: TB-C96LF-SMSM

## Extension Tube



- Latex-free
- Standard 4" (10.7 cm) length
- Optional connectors

Part Number: TB-S4LF-SF

## Inflation System



- Air release valve and inflation bulb
- Durable, latex-free construction

Part Number: AM-INFL-BV

## Replacement Bladders



- Latex-free replacement bladders
- Available in a variety of sizes

Call for Part Number


## Nylons Covers



- Durable nylon construction
- Available in a variety of color coded sizes

Call for Part Number

# Aneroid Part Matrix



Distributed By: