

# Reusable Blood Pressure Cuffs

Two Piece Cuffs – Cover and Bladder



## Features:

- Durable nylon cuff covers
- Comfort fit with no hard edges
- Color coded for easy size identification
- All components are latex-free
- Gauge loop allows for attachment of pocket aneroid
- Arm circumference range guide indicates when larger or smaller cuff should be used
- Full line of connectors allowing compatibility with any NIBP or aneroid
- 3 year warranty

Amico's blood pressure cuffs and components are made with care and quality to ensure they exceed industry standard. Our cuffs are available in all sizes and color coded for quick and easy identification. Our countless connector options are sure to provide the configuration to suit your needs.

Full Range of Sizes



NIBP Compatible



Sphygmomanometer
















# Part Matrix

## Two Tubes



**Reusable Cuffs**

CR2-LFXXX-2QQZZ













<p>Size &amp; Color Options</p> <div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">   <b>07O</b>              Infant – Orange         </div> <div style="text-align: center;">   <b>11N</b>              Adult – Navy Blue         </div> </div> <div style="display: flex; justify-content: space-around; margin-top: 10px;"> <div style="text-align: center;">   <b>09G</b>              Child – Green         </div> <div style="text-align: center;">   <b>12M</b>              Large Adult – Maroon         </div> </div> <div style="display: flex; justify-content: space-around; margin-top: 10px;"> <div style="text-align: center;">   <b>10B</b>              Small Adult – Royal Blue         </div> <div style="text-align: center;">   <b>13W</b>              Thigh – Brown         </div> </div>		<p>Tubes</p>  <b>2</b> Two Tubes	<p>Connector Type (max two)</p> <div style="display: flex; justify-content: space-around; margin-bottom: 10px;"> <div style="text-align: center;">   <b>BV</b>              Inflation Bulb &amp; Valve         </div> <div style="text-align: center;">   <b>BM</b>              Male Bayonet/Rectus         </div> </div> <div style="display: flex; justify-content: space-around; margin-bottom: 10px;"> <div style="text-align: center;">   <b>SM</b>              Male Threaded Screw         </div> <div style="text-align: center;">   <b>QM</b>              Male Quick Connect Subminiature         </div> </div> <div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">   <b>KF</b>              Female Locking Luer         </div> <div style="text-align: center;">   <b>QF</b>              Female Quick Connect Subminiature         </div> </div>
--	--	---	--

## One Tube














**Reusable Cuffs**













CR2-LFXXX-1QQ

<p>Size &amp; Color Options</p> <div style="display: flex; justify-content: space-around;"> <div style="text-align: center;">   <b>07O</b>              Infant – Orange         </div> <div style="text-align: center;">   <b>11N</b>              Adult – Navy Blue         </div> </div> <div style="display: flex; justify-content: space-around; margin-top: 10px;"> <div style="text-align: center;">   <b>09G</b>              Child – Green         </div> <div style="text-align: center;">   <b>12M</b>              Large Adult – Maroon         </div> </div> <div style="display: flex; justify-content: space-around; margin-top: 10px;"> <div style="text-align: center;">   <b>10B</b>              Small Adult – Royal Blue         </div> <div style="text-align: center;">   <b>13W</b>              Thigh – Brown         </div> </div>		<p>Tubes</p>  <b>1</b> One Tube	<p>Connector Type</p> <div style="display: flex; justify-content: space-around; margin-bottom: 10px;"> <div style="text-align: center;">   <b>SM</b>              Male Threaded Screw         </div> <div style="text-align: center;">   <b>QM</b>              Male Quick Connect Subminiature         </div> </div> <div style="display: flex; justify-content: space-around; margin-bottom: 10px;"> <div style="text-align: center;">   <b>KF</b>              Female Locking Luer         </div> <div style="text-align: center;">   <b>QF</b>              Female Quick Connect Subminiature         </div> </div> <div style="text-align: center;">   <b>BM</b>              Male Bayonet/Rectus         </div>
--	--	--	---

# Common NIBP Configurations







 NIBP Side		 Cuff Side		Part Numbers for Two Tube Cuffs
 Male Locking Luer KMKM	 Female Locking Luer KFKF	Infant: CR2-LF070-2KFKF Child: CR2-LF09G-2KFKF Small Adult: CR2-LF10B-2KFKF Adult: CR2-LF11N-2KFKF Large Adult: CR2-LF12M-2KFKF Thigh: CR2-LF13W-2KFKF		
 Female Quick Connect Subminiature QFQF	 Male Quick Connect Subminiature QMQM	Infant: CR2-LF070-2QMQM Child: CR2-LF09G-2QMQM Small Adult: CR2-LF10B-2QMQM Adult: CR2-LF11N-2QMQM Large Adult: CR2-LF12M-2QMQM Thigh: CR2-LF13W-2QMQM		
 Male Quick Connect Subminiature QMQM	 Female Quick Connect Subminiature QFQF	Infant: CR2-LF070-2QFQF Child: CR2-LF09G-2QFQF Small Adult: CR2-LF10B-2QFQF Adult: CR2-LF11N-2QFQF Large Adult: CR2-LF12M-2QFQF Thigh: CR2-LF13W-2QFQF		
 Male/Female Quick Connect Subminiature QMQF	 Female/Male Quick Connect Subminiature QFQM	Infant: CR2-LF070-2QFQM Child: CR2-LF09G-2QFQM Small Adult: CR2-LF10B-2QFQM Adult: CR2-LF11N-2QFQM Large Adult: CR2-LF12M-2QFQM Thigh: CR2-LF13W-2QFQM		
 Female Threaded Screw SFSF	 Male Threaded Screw SMSM	Infant: CR2-LF070-2SMSM Child: CR2-LF09G-2SMSM Small Adult: CR2-LF10B-2SMSM Adult: CR2-LF11N-2SMSM Large Adult: CR2-LF12M-2SMSM Thigh: CR2-LF13W-2SMSM		

 NIBP Side		 Cuff Side		Part Numbers for Single Tube Cuffs
 Female Bayonet/Rectus BF	 Male Bayonet/Rectus BM	Infant: CR2-LF070-1BM Child: CR2-LF09G-1BM Small Adult: CR2-LF10B-1BM Adult: CR2-LF11N-1BM Large Adult: CR2-LF12M-1BM Thigh: CR2-LF13W-1BM		
 Male Locking Luer KM	 Female Locking Luer KF	Infant: CR2-LF070-1KF Child: CR2-LF09G-1KF Small Adult: CR2-LF10B-1KF Adult: CR2-LF11N-1KF Large Adult: CR2-LF12M-1KF Thigh: CR2-LF13W-1KF		
 Female Quick Connect Subminiature QF	 Male Quick Connect Subminiature QM	Infant: CR2-LF070-1QM Child: CR2-LF09G-1QM Small Adult: CR2-LF10B-1QM Adult: CR2-LF11N-1QM Large Adult: CR2-LF12M-1QM Thigh: CR2-LF13W-1QM		
 Male Quick Connect Subminiature QM	 Female Quick Connect Subminiature QF	Infant: CR2-LF070-1QF Child: CR2-LF09G-1QF Small Adult: CR2-LF10B-1QF Adult: CR2-LF11N-1QF Large Adult: CR2-LF12M-1QF Thigh: CR2-LF13W-1QF		
 Female Threaded Screw SF	 Male Threaded Screw SM	Infant: CR2-LF070-1SM Child: CR2-LF09G-1SM Small Adult: CR2-LF10B-1SM Adult: CR2-LF11N-1SM Large Adult: CR2-LF12M-1SM Thigh: CR2-LF13W-1SM		

# Blood Pressure Cuff Components



		Limb Size	Cuff Cover	Two Tube Bladder	One Tube Bladder
Infant	 Orange	9 - 13 cm (3.5 - 5 in)	CR2-CVR07O	CR2-BLDR07R-2	CR2-BLDR07R-1
Child	 Green	12 - 20 cm (4.7 - 7.9 in)	CR2-CVR09G	CR2-BLDR09R-2	CR2-BLDR09R-1
Small Adult	 Royal Blue	19 - 27 cm (7.5 - 10.6 in)	CR2-CVR10B	CR2-BLDR10R-2	CR2-BLDR10R-1
Adult	 Navy Blue	26 - 36 cm (10.2 - 14.2 in)	CR2-CVR11N	CR2-BLDR11R-2	CR2-BLDR11R-1
Large Adult	 Maroon	34 - 47 cm (13.4 - 18.5 in)	CR2-CVR12M	CR2-BLDR12R-2	CR2-BLDR12R-1
Thigh	 Brown	45 - 57 cm (17.7 - 22.4 in)	CR2-CVR13W	CR2-BLDR13R-2	CR2-BLDR13R-1

## Wall Aneroid



- Accurate and easy to read
- Independent cuff basket allows proper placement for ease of access
- Variety of rail adapters

Part Number: AM-AB-LF2118

## Pocket Aneroid



- Compact size
- Sturdy construction
- Sold with or without inflation system

## Stethoscopes



- Stainless steel double-head
- Non-chill rim
- Adult and pediatric available
- Latex-free

Call for Part Number

## Extension Tube



- Latex-free tubing
- Optional connectors

Please Call for Part Number

## Inflation Bulb



- Air release valve and inflation bulb
- Durable, latex-free construction

Part Number: AM-INFL-BV

## Connectors



- Large variety of connectors and tubing

## Baskets



- Various basket sizes available to contain blood pressure cuffs

## Distributed By:

