

Diagnostic Equipment Catalog



Amico



Dedicated to Patient
and Practitioner Experience



Contents

Instrument Heads

Ophthalmoscope	5
Coaxial	5
Coaxial Precision	5
Otoscope	6
Fiber Optic Otoscope	6
Operating Otoscope	6
Streak Retinoscope	7
Illuminators	7
Nasal Illuminator	7
Bi-valve Nasal Speculum	7
Larynx	7
Fiber-Optic Bent	7
Consumables	8
Bulbs	8
Specula	9

Diagnostic Station

Direct Wall	11
Direct Wall Common Configurations	12
Wallboard Mounted	13
Wallboard Mounted Common Configurations	14
Wallboard Part Numbers	14
Rail Mounted	15
Rail Mounted Common Configurations	16-17
Ordering Guide	17
Rail Options	18
Rail Accessories	18
Rail Mounted with Integrated Storage	19
Rail Mounted with Integrated Configurations	20
Roll Stand	21
Roll Stand Mounted Common Configurations	22-23

Aneroid/Sphygmomanometer

Wall Aneroid/Sphygmomanometer	25
Mobile Aneroid	26
Rail Mount Aneroid	27
Aneroid Accessories Inflation System Parts	28

Blood Pressure/Cuffs

Two-piece Reusable Cuffs	31
One-tube Cuff Configurations	32
Two-tube Cuff Configurations	32
Blood Pressure Cuff Connectors	33



01

Instrument Heads

Ophthalmoscope



The Amico Ophthalmoscope provides superior optics and is available in Coaxial and Coaxial Precision models.

- | | | | |
|--|-------------------|--|---------------------|
| | Small Aperture | | Cobalt Blue Filter |
| | Medium Aperture | | Polarization Filter |
| | Large Aperture | | Red-free Filter |
| | Slit Aperture | | Unfiltered |
| | Fixation Aperture | | |

Coaxial



- 28 adjustable diopter lenses (+40 to -25)
- 18 unique aperture and filter combinations
- Sealed design keeps dust out
- Anti-theft locking feature
- Halogen and LED available
- Three (3) year warranty

Part Number	Bulb
DH-W35-CH-L	Halogen
DH-W35-CL-L	LED

LIFETIME Warranty on LED Bulb included with this product, when you order the LED model.

Coaxial Precision



- 68 single-diopter steps (+38 to -30) for precise focus
- 18 unique aperture and filter combinations
- Sealed design keeps dust out
- Anti-theft locking feature
- Halogen and LED available
- Three (3) year warranty

Part Number	Bulb
DH-W35-PH-L	Halogen
DH-W35-PL-L	LED

LIFETIME Warranty on LED Bulb included with this product, when you order the LED model.

Otoscope

Fiber Optic Otoscope



The Amico Otoscope offers an ultra-clear field of view, with swivel lens and secure lock included standard.

- Fiber optics project light without obstruction
- Swivelling lens allows instrument access
- Compatible with insufflation bulb
- Anti-theft locking feature
- Halogen and LED available
- Three (3) year warranty

Part Number	Bulb
DH-W35-FH-L	Halogen
DH-W35-FL-L	LED

Insufflation Bulb



- Otoscope insufflation bulb and tube
- Compatible with fiber optic otoscope

Part Number
DA-OS-IB

Operating Otoscope



- Fiber optic illumination
- Design allows unhindered access to ear
- Swivelling lens and adjustable speculum holder
- Three (3) year warranty

Part Number	Bulb
DH-W35-OH	Halogen
DH-W35-OL	LED

Streak Retinoscope

Retinoscope



- Single control sleeve for vergence and rotation
- Superior streak quality enables easier refractions
- Sealed design keeps internal components dust free
- Anti-theft locking feature
- Three (3) year warranty

Part Number	Bulb
DH-W35-RH-L	Halogen
DH-W35-RL-L	LED

Precision is always kept in mind with all of Amico Diagnostic products.

Illuminators

Nasal



Nasal Illuminator, 3.5 V. Head Only.

Part Number	Bulb
DH-W35-NH	Halogen
DH-W35-NL	LED

- Three (3) year warranty

Bi-valve Nasal Speculum



Bi-valve Nasal Speculum Illuminator, 3.5 V. Head Only.

Part Number	Bulb
DH-W35-VH	Halogen
DH-W35-VL	LED

- Three (3) year warranty

Larynx



Larynx Illuminator, 3.5 V. Head Only.

Part Number	Bulb
DH-W35-LH	Halogen
DH-W35-LL	LED

- Three (3) year warranty

Fiber-Optic Bent



Bent Illuminator, 3.5 V. Head Only.

Part Number	Bulb
DH-W35-BH	Halogen
DH-W35-BL	LED

- Three (3) year warranty

Consumables

Bulbs

Amico offers halogen and LED replacement bulbs to ensure you have the brightness and bulb-life to meet your needs.

- LED bulbs use approximately 85.7% less power and have 1,000x the lifespan of a halogen bulb. Save money on consumables with LEDs
- LED Bulbs are brighter than halogen bulbs

Lifetime warranty on LED bulbs when you order an LED Diagnostic Station.



View of LED Bulb scope

LED



Part Number	Application and Packaging
DC-3100-L01	LED bulb 3100 series for Fiber optic Otoscope and Illuminators, pack of 1
DC-4900-L01	LED bulb 4900 series for Coaxial Ophthalmoscope, pack of 1
DC-1100-L01	LED bulb 1100 series for Retinoscope, pack of 1

Halogen



Part Number	Application and Packaging
DC-3100-H10	Halogen bulb 3100 series for Fiber optic Otoscope, pack of 10
DC-4900-H10	Halogen bulb 4900 series for Coaxial Ophthalmoscope, pack of 10
DC-1100-H01	Halogen bulb 1100 series for Retinoscope, pack of 1

Specula

Fiber Optic Otoscope Speculas



Size	Disposable
● 2.5 mm	DC-SSD-1K-25M, Pack of 1,000
● 3.0 mm	DC-SSD-1K-30M, Pack of 1,000
● 4.0 mm	DC-SSD-1K-40M, Pack of 1,000
● 5.0 mm	DC-SSD-1K-50M, Pack of 1,000

Operating Otoscope Speculas



Size	Disposable
● 2.5 mm	DO-OPD-1K-25M, Pack of 1,000
● 3.0 mm	DO-OPD-1K-30M, Pack of 1,000
● 4.0 mm	DO-OPD-1K-40M, Pack of 1,000
● 5.0 mm	DO-OPD-1K-50M, Pack of 1,000





02
Diagnostic Station

Direct Wall



Smart practices and facilities are always searching for new ways to eliminate headaches, cut costs, and optimize processes. Developed in collaboration with doctors and nurses, Amico's diagnostic equipment line maximizes value, provides peace of mind and enhances practitioner workflow.

A range of mounting options optimize workflow and enhance practitioner and patient experience. Mount directly to the wall, or on a wallboard to enhance aesthetics and simplify mounting. Roll Stand mounting enables quick relocation from room to room in a busy office or handling patient overflow in a hectic emergency department. Rails facilitate easy mounting, quick reorganization of equipment, and incorporates even more storage for the ultimate exam workflow solution!

Features



Standard Instrument Heads Include Anti-theft Locks



Integrated Cable Management



USB Outlets for Auxiliary Devices



Locking Thermometer














Integrated Dispenser Storage



Eco-friendly Low Power Consumption

Direct Wall Common Configurations

		Coaxial Ophthalmoscope	Fiber Optic Otoscope	Specula Dispenser	Aneroid with Adult Cuff and Basket	Thermometry	
						Oral/Axillary	Infrared Tympanic
							
Configuration	Part Number						
	DS-UDM-CHFHH-XXX	Halogen	Halogen				
	DS-UDM-CHFL-XXX	Halogen	LED				
	DS-UDM-CLFL-XXX	LED	LED				
	DS-UDM-CHFHH-DXX	Halogen	Halogen				
	DS-UDM-CHFL-DXX	Halogen	LED	◆			
	DS-UDM-CLFL-DXX	LED	LED				
	DS-UDM-CHFHH-DAX	Halogen	Halogen				
	DS-UDM-CHFL-DAX	Halogen	LED	◆	◆		
	DS-UDM-CLFL-DAX	LED	LED				
	DS-UDM-CHFHH-DAO	Halogen	Halogen				
	DS-UDM-CHFL-DAO	Halogen	LED	◆	◆	◆	
	DS-UDM-CLFL-DAO	LED	LED				
	DS-UDM-CHFHH-DAI	Halogen	Halogen				
	DS-UDM-CHFL-DAI	Halogen	LED	◆	◆		◆
	DS-UDM-CLFL-DAI	LED	LED				

Wallboard Mounted



Avoid installation pitfalls by using Amico Wallboards to securely fasten diagnostic equipment to your wall.

- Decorative thermofoil over wood, for easy cleaning, durability and better esthetics
- Power cord concealed behind board
- Custom shapes for specialty areas
- Pre-drilled wallboard allows for simple installation of diagnostic station
- Wide variety of colors and wood laminates to meet all design needs
- Integrated cable management

Color Options




Light Grey (Standard)



Candlelight (Optional)



Medium Oak (Optional)



Hardrock Maple (Optional)



Concealed Cable Management

Specialty Wallboards



Blue Whale












Green Dinosaur





Yellow Bus


Wallboard Mounted Common Configurations

		Coaxial Ophthalmoscope	Fiber Optic Otoscope	Specula Dispenser	Aneroid with Adult Cuff and Basket	Thermometry	
						Oral/Axillary	Infrared Tympanic
							
Configuration	Part Number						
	DS-ULG-CHFh-DAX	Halogen	Halogen				
	DS-ULG-CHFL-DAX	Halogen	LED	◆	◆		
	DS-ULG-CLFL-DAX	LED	LED				
	DS-ULG-CHFh-DAO	Halogen	Halogen				
	DS-ULG-CHFL-DAO	Halogen	LED	◆	◆	◆	
	DS-ULG-CLFL-DAO	LED	LED				
	DS-ULG-CHFh-DAI	Halogen	Halogen				
	DS-ULG-CHFL-DAI	Halogen	LED	◆	◆		◆
	DS-ULG-CLFL-DAI	LED	LED				

Wallboard Part Numbers

	Part Number	Color	Size
	DS-WLLBRD-LG03	LG - Light Grey	For Diagnostic Station with Dispenser and Aneroid
	DS-WLLBRD-LG04	LG - Light Grey	For Diagnostic Station with Dispenser Aneroid and Thermometry

Rail Mounted



- Easily adjust and relocate equipment to suit the needs of your practitioners
- Equipment can be locked in place
- Simple transfer from direct rail to headwall or casework rail
- Improve patient-centered care and HCAHPS Scores
- Facilitates terminal cleaning

Features and Configurations



Vertical, Horizontal, and Fairfield Rail and Mounts















Easily Relocate Your Equipment



Custom Configuration with Storage

Rail Mounted Common Configurations

			Amico Horizontal Rail Mount	Coaxial Ophthalmoscope	Fiber Optic Otoscope	Specula Dispenser	Aneroid with Adult Cuff and Basket	Aneroid on L-Bracket with Basket	Thermometry	
									Oral/ Axillary	Infrared Tympanic
										
*Configuration	Part Number									
	DS-UA1-CHFh-XXX	◆	Halogen	Halogen						
	DS-UA1-CHFL-XXX		Halogen	LED						
	DS-UA1-CLFL-XXX		LED	LED						
	DS-UA1-CHFh-DXX	◆	Halogen	Halogen						
	DS-UA1-CHFL-DXX		Halogen	LED	◆					
	DS-UA1-CLFL-DXX		LED	LED						
	DS-UA1-CHFh-DAX	◆	Halogen	Halogen						
	DS-UA1-CHFL-DAX		Halogen	LED	◆	◆				
	DS-UA1-CLFL-DAX		LED	LED						
	DS-UA1-CHFh-DLX	◆	Halogen	Halogen						
	DS-UA1-CHFL-DLX		Halogen	LED	◆		◆			
	DS-UA1-CLFL-DLX		LED	LED						
	DS-UA1-CHFh-DAO	◆	Halogen	Halogen						
	DS-UA1-CHFL-DAO		Halogen	LED	◆	◆		◆		
	DS-UA1-CLFL-DAO		LED	LED						




* Rail not included. For rail information, see page 18.

Rail Mounted Common Configuration

Amico Horizontal Rail Mount	Coaxial Ophthalmoscope	Fiber Optic Otoscope	Specula Dispenser	Aneroid with Adult Cuff and Basket	Aneroid on L-Bracket with Basket	Thermometry	
						Oral/ Axillary	Infrared Tympanic
							

*Configuration

Part Number

	DS-UA1-CHFHDLO		Halogen	Halogen					
	DS-UA1-CHFLDLO	◆	Halogen	LED	◆		◆	◆	
	DS-UA1-CLFLDLO		LED	LED					
	DS-UA1-CHFHDAI		Halogen	Halogen					
	DS-UA1-CHFLDAI	◆	Halogen	LED	◆	◆			◆
	DS-UA1-CLFLDAI		LED	LED					
	DS-UA1-CHFHDLI		Halogen	Halogen					
	DS-UA1-CHFLDLI	◆	Halogen	LED	◆		◆		◆
	DS-UA1-CLFLDLI		LED	LED					

* Rail not included. For rail information, see page 18.

Did you know you can combine diagnostics and storage for a complete workflow solution? See page 19.

Ordering Guide

Do you need mounts for a different type of rail? Use the ordering guide below to determine your part number or ask your sales representative for assistance.

Other rail mount options available.
Add your code in the ordering guide.

For example: DS-UF1-CHFHDAX

Diagnostic Station, Fairfield rail mounts, Halogen Bulb, Specula Dispenser, Aneroid

DS-UA1-CHFHDAX




A1 Horizontal Rail Mounts

V1 Vertical Rail Mounts

F1 Fairfield Rail Mounts

Rail Options


Amico Horizontal Rail System



Part Number	Length
ARS-024	24"
ARS-036	36"
ARS-048	48"

- Integrated mounting plate
- Ease of maintenance and infection control
- Available in a wide variety of lengths
- Laminate insert allows for color customization


Amico Vertical Rail System



Part Number	Length
VRS-024	24"
VRS-036	36"
VRS-048	48"

- Mounts flush to any surface
- Can be used in conjunction with horizontal rails for a more efficient use of space
- Ideal for inlaying in cabinets/casework


Light Duty Fairfield Rail System



Part Number	Length
LRS-024	24"
LRS-036	36"
LRS-048	48"

- Slim profile ideal for adding to existing equipment
- Laminate insert allows for color customization
- Available in a wide variety of lengths

Fairfield Rail System



Part Number	Length
FRS-024	24"
FRS-036	36"
FRS-048	48"

- Can be fitted directly to the wall surface or paired with mounting hardware
- Laminate insert allows for color customization
- Available in a wide variety of lengths

Rail Accessories

Basket 5" x 5" x 4"

(12.7 cm x 12.7 cm x 10.2 cm)



Small storage baskets are great for BP cuffs and miscellaneous odds and ends. Designed to sit below the horizontal rail or in the center for vertical systems.

Part Number	Type
050504-A-BB	Amico Horizontal Rail
050504-V-BB	Amico Vertical Rail
050504-F-BB	Fairfield Rail

Basket 9" x 9" x 10"

(22.9 cm x 22.9 cm x 25.4 cm)



Medium is designed to accommodate small objects for nurses and doctors convenience. Designed to sit below the horizontal rail or in the center for vertical rail systems.

Part Number	Type
090910-A-BB	Amico Horizontal Rail
090910-V-BB	Amico Vertical Rail
090910-F-BB	Fairfield Rail

Basket 10" x 6" x 9"

(25.4 cm x 15.2 cm x 22.9 cm)



This basket works well with small items in patient rooms can also work as a convenient waste receptacle.

Part Number	Type
100609-A-BB	Amico Horizontal Rail
100609-V-BB	Amico Vertical Rail
100609-F-BB	Fairfield Rail

Glove/Tissue Box Holder



Versatile storage for glove boxes or tissue boxes, can be mounted onto various style rails or direct to wall.

Part Number	Type
GLBOXHLD-1-A	Amico Horizontal Rail
GLBOXHLD-1-V	Amico Vertical Rail
GLBOXHLD-1-F	Fairfield Rail

Rail Mounted with Integrated Storage



Have something you'd like to mount not shown on this page?
Ask your sales rep about custom mounting!

- Faster workflow
- Better organization
- Industry-best warranty
- Superb storage
- Customizable
- Extendable
- More patient face time

Key Storage Features



Wire Basket

Great for storing extra cuffs and medium/large supplies.



Glove/Tissue Box Holder

Fits most sizes of glove and tissue boxes. Spring loaded to accommodate shallow boxes.



Stethoscope Hook

Great for storing your stethoscope next to your aneroid.



Clear Storage Container

Perfect for your swabs, tongue depressors, cotton balls, and more!



Sharps Container Holder















Universal Sharps holder with adjustable strap fits common sizes of sharp bins.



4' Rail

4' Fairfield Rail for quick mounting of your diagnostic station and storage.

Rail Mounted with Integrated Configuration

		DIAGNOSTICS					STORAGE				
		4' Fairfield Rail	Coaxial Ophthalmoscope	Fiber Optic Otoscope	Specula Dispenser	Aneroid with Adult Cuff and Basket	Oral/Axillary Thermometer	Wire Basket	Clear Storage Container	Glove/Tissue Box Holder	Stethoscope Hook
Configuration	Part Number										
	DS-UF3-CHFH-DAX		Halogen	Halogen							
	DS-UF3-CHFL-DAX	◆	Halogen	LED	◆	◆			◆ ◆	◆	◆
	DS-UF3-CLFL-DAX		LED	LED							
	DS-UF3-CHFH-DAO		Halogen	Halogen							
	DS-UF3-CHFL-DAO	◆	Halogen	LED	◆	◆	◆		◆ ◆	◆	◆
	DS-UF3-CLFL-DAO		LED	LED							
	DS-UF4-CHFH-DAX		Halogen	Halogen							
	DS-UF4-CHFL-DAX	◆	Halogen	LED	◆	◆		◆ ◆		◆	◆
	DS-UF4-CLFL-DAX		LED	LED							
	DS-UF4-CHFH-DAO		Halogen	Halogen							
	DS-UF4-CHFL-DAO	◆	Halogen	LED	◆	◆	◆	◆ ◆		◆	◆
	DS-UF4-CLFL-DAO		LED	LED							

Roll Stand

Modular Design Adds Flexibility

Mounts Compatible with Amico Horizontal Rails

Mount your diagnostic system on a Roll Stand for complete flexibility. Easily transport your station from room to room to suit the dynamic needs of your facility.

- Modular design allows devices to be easily added and removed
- Device mounts are compatible with Amico horizontal headwalls
- 10 ft (3 m) hospital grade power cord
- Extra large basket for storage
- All casters are locking
- Stable 5-leg recessed Base
- Available with thermometry and aneroid
- 10 ft (3 m) power cord











Large Storage Basket

Ultra-stable Recessed Base

All Locking Wheels

Perfect for patient overflow!

Roll Stand Mounted Common Configurations

Configuration	Part Number						Thermometry	
							Oral/Axillary	Infrared Tympanic
	DS-UR1-CHFHDXX	◆	Halogen	Halogen	◆			
	DS-UR1-CHFLDXX		Halogen	LED				
	DS-UR1-CLFLDXX		LED	LED				
	DS-UR1-CHFHDAX	◆	Halogen	Halogen	◆	◆		
	DS-UR1-CHFLDAX		Halogen	LED				
	DS-UR1-CLFLDAX		LED	LED				
	DS-UR1-CHFHD AO	◆	Halogen	Halogen	◆	◆	◆	
	DS-UR1-CHFLD AO		Halogen	LED				
	DS-UR1-CLFLD AO		LED	LED				
	DS-UR1-CHFHD AI	◆	Halogen	Halogen	◆	◆		◆
	DS-UR1-CHFLD AI		Halogen	LED				
	DS-UR1-CLFLD AI		LED	LED				
	DS-UR1-CHFHD XO	◆	Halogen	Halogen	◆		◆	
	DS-UR1-CHFLD XO		Halogen	LED				
	DS-UR1-CLFLD XO		LED	LED				

Roll Stand Mounted Common Configurations

Configuration	Part Number	Roll Stand	Coaxial Ophthalmoscope	Fiber Optic Otoscope	Specula Dispenser	Aneroid and Cuff Basket	Thermometry	
							Oral/Axillary	Infrared Tympanic
								



DS-UR1-CHFHD-DXI	◆	Halogen	Halogen	◆			◆
DS-UR1-CHFL-DXI		Halogen	LED				
DS-UR1-CLFL-DXI		LED	LED				



Model Shown: DS-UR1-CHFL-DAI



Aneroid/
Sphygmomanometer

03

Wall Aneroid/Sphygmomanometer

6" (15.2 cm) Luminescent Dial



Large Independent Cuff Basket

Latex-free Inflation System



Amico's BP Aneroid Sphygmomanometers provide an accurate pressure measurements. The clarity of the dial makes it easy to read, and features a scratch and shatter resistant lens. Available with any cuff size or fittings to match any blood pressure cuff.

- Luminescent dial is readable in low light
- Large 6" (15.2 cm) dial
- Recessed dial face reduces parallax errors
- Lens minimizes reflection to diminish distracting mirror effect
- Thermally treated diaphragm maintains calibration
- 110° gauge swivel for easy viewing
- Available with latex-free blood pressure cuff and 8 ft (2.4 m) coiled tubing inflation system
- No pin stop movement
- Mercury-free
- Lifetime calibration warranty

Part Number

AM-DB-LF2118

See page 31 for extra Cuff sizes.

Huge Value!

Mobile Aneroid



The Mobile Aneroid enables quick and easy movement from room to room.

- Luminescent dial is readable in low light
- Accurate to within 3 mmHg
- Large 6" (15.2 cm) dial
- Recessed dial face reduces parallax errors
- Lens minimizes reflection to diminish distracting mirror effect
- 110° gauge swivel for easy viewing
- Mercury free and latex-free
- Adjustable height from 34" to 48" (86.4 cm to 121.9 cm)
- 18" (45.7 cm) diameter weighted base for added stability
- Spring-loaded post helps prevent equipment damage
- Aneroid face angled upwards at 20° for ease of viewing

Part Number

AM-RX-LF2118

See page 31 for more information on Cuffs.

Rail Mount Aneroid



Mount your aneroid, cuff basket, and other diagnostic products using Amico's rail system for simplified installation, enhanced workflow, and easy modification of equipment and storage layouts.

- Reduces congestion on the equipment rails
- Available with a variety of rail adapters

See page 20 for more information on rail-mounted configurations.

See page 31 for information about ordering extra Cuffs.

Rail Mounts



Model Shown: AM-AB-LF2118

Model Shown: AM-VB-LF2118

- Easy to move accessory adapters
- Independent cuff basket allows proper placement for ease of access
- Available with a variety of rail adapters

Part Number	Rail Mount
AM-AB-LF2118	Amico Horizontal
AM-VB-LF2118	Amico Vertical
AM-FB-LF2118	Fairfield

L - Bracket Mount



Model Shown: AM-AL-LF2118

- Aneroid and basket mounted together on a single bracket
- Reduces congestion on the equipment rails
- Available with a variety of rail adapters

Part Number	Rail Mount
AM-AL-LF2118	Amico Horizontal
AM-VL-LF2118	Amico Vertical
AM-FL-LF2118	Fairfield

Aneroid Accessories and Inflation System Parts

Coil Tubing



Latex-free coil tubing.
Expands to 8 ft (2.4 m).

Part Number
TB-C96LF

Coil Tubing & Extension



8 ft latex-free coil tubing with
4" extension and male screw connector.

Part Number
AM-COILASM8-SM

Inflation System



Inflation bulb with air release and end valve.
Durable, latex-free construction.

Part Number
AM-INFL-BV

Lens Cover



Replacement lens for wall aneroid.

Part Number
AM-LENS

Aneroid Gauge



Aneroid gauge only, 300 mmHG

Part Number	Mount
AM-A-SPHYG	Amico Horizontal Rail
AM-V-SPHYG	Amico Vertical Rail
AM-F-SPHYG	Fairfield Rail
AM-D-SPHYG	Wall

Cuff Basket



Low profile cuff storage.

Part Number	Mount
AM-A-CUFBSK	Amico Horizontal Rail
AM-V-CUFBSK	Amico Vertical Rail
AM-F-CUFBSK	Fairfield
AM-D-CUFBSK	Wall





Ask about our other products such as the **Hummingbird Workstation on Wheels** with **integrated Diagnostic Station!**



04

Blood Pressure Cuffs

Two-piece Reusable Cuffs









Choose your Cuffs for your
Amico Aneroid!

Amico's blood pressure cuffs and components are made with care and quality to ensure they exceed industry standard. Our cuffs are available in all sizes and color coded for quick and easy identification.

- Durable nylon cuff covers
- Comfort fit with no hard edges
- Color coded for easy size identification
- All components are latex-free
- Arm circumference range guide indicates when larger or smaller cuff should be used
- Full line of connectors allowing compatibility with any NIBP or aneroid
- 3 year warranty



Two-piece cuff with
bulb and connector.

Size	Color	Circumference	Part Number	Description
Infant	 Orange	9 - 13 cm (3.5 - 5 in)	CR2-LF07O-2BVSF	Two-piece Reusable Blood Pressure Cuff. Inflation Bulb and Female Threaded Screw Connector.
Child	 Green	12 - 20 cm (4.7 - 7.9 in)	CR2-LF09G-2BVSF	Two-piece Reusable Blood Pressure Cuff. Inflation Bulb and Female Threaded Screw Connector.
Small Adult	 Royal Blue	19 - 27 cm (7.5 - 10.6 in)	CR2-LF10B-2BVSF	Two-piece Reusable Blood Pressure Cuff. Inflation Bulb and Female Threaded Screw Connector.
Adult	 Navy Blue	26 - 36 cm (10.2 - 14.2 in)	CR2-LF11N-2BVSF	Two-piece Reusable Blood Pressure Cuff. Inflation Bulb and Female Threaded Screw Connector.
Large Adult	 Maroon	34 - 47 cm (13.4 - 18.5 in)	CR2-LF12M-2BVSF	Two-piece Reusable Blood Pressure Cuff. Inflation Bulb and Female Threaded Screw Connector.
Thigh	 Brown	45 - 57 cm (17.7 - 22.4 in)	CR2-LF13W-2BVSF	Two-piece Reusable Blood Pressure Cuff. Inflation Bulb and Female Threaded Screw Connector.

One-tube Cuff Configurations



Cuff with one-tube,
Male Bayonet Connector.

CR2-LF07O-1BM	Infant
CR2-LF09G-1BM	Child
CR2-LF10B-1BM	Small Adult
CR2-LF11N-1BM	Adult
CR2-LF12M-1BM	Large
CR2-LF13W-1BM	Thigh

Cuff with one-tube,
No Connector.

CR2-LF07O-1	Infant
CR2-LF09G-1	Child
CR2-LF10B-1	Small Adult
CR2-LF11N-1	Adult
CR2-LF12M-1	Large
CR2-LF13W-1	Thigh

Cuff with one-tube,
Female Locking Luer Connector.

CR2-LF07O-1KF	Infant
CR2-LF09G-1KF	Child
CR2-LF10B-1KF	Small Adult
CR2-LF11N-1KF	Adult
CR2-LF12M-1KF	Large
CR2-LF13W-1KF	Thigh

Cuff with one-tube,
Male Subminiature Connector.

CR2-LF07O-1QM	Infant
CR2-LF09G-1QM	Child
CR2-LF10B-1QM	Small Adult
CR2-LF11N-1QM	Adult
CR2-LF12M-1QM	Large
CR2-LF13W-1QM	Thigh

Cuff with one-tube,
Female Subminiature Connector.

CR2-LF07O-1QF	Infant
CR2-LF09G-1QF	Child
CR2-LF10B-1QF	Small Adult
CR2-LF11N-1QF	Adult
CR2-LF12M-1QF	Large
CR2-LF13W-1QF	Thigh

Cuff with one-tube,
Male Screw Connector.

CR2-LF07O-1SM	Infant
CR2-LF09G-1SM	Child
CR2-LF10B-1SM	Small Adult
CR2-LF11N-1SM	Adult
CR2-LF12M-1SM	Large
CR2-LF13W-1SM	Thigh



Cuff with two-tubes,
No Connectors.

CR2-LF07O-2	Infant
CR2-LF09G-2	Child
CR2-LF10B-2	Small Adult
CR2-LF11N-2	Adult
CR2-LF12M-2	Large Adult
CR2-LF12M-2	Thigh

Cuff with two-tubes,
Inflation Bulb & Male Screw Connector.

CR2-LF07O-SMSM	Infant
CR2-LF09G-SMSM	Child
CR2-LF10B-SMSM	Small Adult
CR2-LF11N-SMSM	Adult
CR2-LF12M-SMSM	Large
CR2-LF13W-SMSM	Thigh

Cuff with two-tubes,
Female Locking Luer Connectors.

CR2-LF07O-2KFKF	Infant
CR2-LF09G-2KFKF	Child
CR2-LF10B-2KFKF	Small Adult
CR2-LF11N-2KFKF	Adult
CR2-LF12M-2KFKF	Large
CR2-LF13W-2KFKF	Thigh

Cuff with two-tubes,
Mating Subminiature Connectors.

CR2-LF07O-2QFQM	Infant
CR2-LF09G-2QFQM	Child
CR2-LF10B-2QFQM	Small Adult
CR2-LF11N-2QFQM	Adult
CR2-LF12M-2QFQM	Large
CR2-LF13W-2QFQM	Thigh

Cuff with two-tubes,
Female Subminiature Connectors.
















CR2-LF07O-2QFQF	Infant
CR2-LF09G-2QFQF	Child
CR2-LF10B-2QFQF	Small Adult
CR2-LF11N-2QFQF	Adult
CR2-LF12M-2QFQF	Large
CR2-LF13W-2QFQF	Thigh

Cuff with two-tubes, Inflation Bulb
& Female Locking Luer Connector.

CR2-LF07O-2BVKF	Infant
CR2-LF09G-2BVKF	Child
CR2-LF10B-2BVKF	Small Adult
CR2-LF11N-2BVKF	Adult
CR2-LF12M-2BVKF	Large
CR2-LF13W-2BVKF	Thigh

Two-tube Cuff Configurations

Blood Pressure Cuff Connectors

TC-SFP-532-10	Female Screw with Plastic Connector 5/32" ID Tubing 10/Pkg.		TC-KFP-332-10	Female Locking Luer with Plastic Connector 3/32" ID Hose 10/Pkg.	
TC-SMP-532-10	Male Screw with Plastic Connector 5/32" ID Tubing 10/Pkg.		TC-KMP-532-10	Male Locking Luer Connector 5/32" ID Hose 10/Pkg.	
TC-QFP-532-10	Female Subminiature with Plastic Connector 5/32" ID Hose 10/Pkg.		TC-LFP-532-10	Female Slip Luer with Plastic Connector 5/32" ID Hose 10/Pkg.	
TC-QMP-532-10	Male Subminiature with Plastic Connector 5/32" ID Hose 10/Pkg.		TC-LMP-532-10	Male Slip Luer with Plastic Connector 5/32" ID Hose 10/Pkg.	
TC-BFM-532	Female Bayonet with Metal Connector for 5/32" ID Tubing.		TC-LMP-332-10	Male Slip Luer with Plastic Connector 3/32" ID Hose 10/Pkg.	
TC-BMM-532	Male Bayonet with Metal Connector for 5/32" ID Tubing.		TC-HBP-532-10	Double Ended Hose Barb Connector for 5/32" ID Hose 10/Pkg.	
TC-BMP-532-10	Male Bayonet with Plastic Connector for 5/32" ID Tubing 10/Pkg.		TC-HBP-332-10	Double Ended Hose Barb Connector for 3/32" ID Hose 10/Pkg.	
TC-KFP-532-10	Locking Luer with Plastic Connector 5/32" ID Hose 10/Pkg.				



*A*mico

Right at your finger tips.
Patient exams made easy.



Instrument Heads



Diagnostic Station



Aneroid/ Sphygmomanometer



Blood Pressure Cuffs





Amico Accessories Inc. | 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada
Toll Free Tel: 1.877.264.2697 | Tel: 905.763.7778 | Fax: 905.763.8587
Email: accessories@amico.com | www.amico.com

Represented By: