

Diagnostic Station

Modular System



Part Number: DS-UDM-CHFH-DXX

Features

- Modular design allows for easy customization
- Ergonomic handle with automatic ON/OFF function
- Low power, long life LED options
- High efficiency power supply reduces energy costs
- Smart Bulb System with timer control extends bulb life and saves energy
- Space saving with smallest equipment footprint available
- Rheostat on handle allows easy control of light intensity
- Handle's spiral cord expands to over 10 ft (3 m)
- Easy to clean and disinfect
- Comes standard with anti-theft locking instrument heads
- Latex and mercury free aneroid and cuff
- Multiple mounting possibilities to suit specific needs
- Compatible with majority of manufacturers' instrument heads

Amico Diagnostic Incorporated provides quality manufactured medical equipment. The diagnostic station with Smart Bulb system saves on energy, consumables and space. Our countless mounting options will provide a solution for any workflow.



Part Number: DS-UA1-CHFH-DAO

Mounting and Module Options

Mounting Options


Wall Board



- Pre-drilled wall board allows for simple installation of diagnostic station
- Wide variety of colors and wood laminates to meet all design needs
- Color can be matched to Amico Headwalls
- Cable management behind board

Part Number: DS-ULG-CHFH-DAO

Specialty Wall Boards



- Helps to provide non-clinical ambiance
- Pre-drilled board with cable management behind board
- Customized shapes available

Call for part number


Roll Stand



- Large 5 leg base for increased stability
- Large basket allows for easy storage of accessories
- Cord wrap for proper management of power cord
- Option to accommodate thermometer and aneroid

Part Number: DS-UR1-CHFH-DXX

Rail Mount



- More flexible positioning over traditional wall board
- Equipment can be locked down or moved easily
- Can be mounted on any headwall, flatwall or direct rail option

Call for part number

Module Options

Base Station



- Automatic ON/OFF independent of handle cradle
- High efficiency power supply allows simultaneous instrument operation
- Integrated cord wrap eliminates excess wiring
- Timer shut-off preserves bulb life



Part Number: DS-UDM-2

Specula Dispenser



- Dispenses up to 4 specula sizes
- Compartment for deep or shallow storage
- Clear housing allows for quick inventory management
- Compatible with leading brand of specula

Part Number: DS-DCD-SPDISP

Wall Aneroid



- Accurate within 3 mmHg
- 110° swivel adjustment
- Large 6" (15 cm) recessed dial face for increased visibility
- Latex and mercury free
- Lifetime calibration warranty

Part Number: AM-DB-LF2118

Electronic Thermometry



- Predictive or infrared options available
- Fast read response in under 4 seconds
- 60 second pulse timer option
- Hygienic single use probe covers allow touchless loading and ejecting

Part Number: TT-D-IR

Instrument Head Options

Fiber Optic Otoscope



- Fiber optics project light without obstruction
- Swiveling lens allows instrument access
- Compatible with insufflation bulb
- Anti-theft locking feature
- Halogen and LED available

Part Number: DH-W35-FL-L: LED
Part Number: DH-W35-FH-L: Halogen

Operating Otoscope



- Fiber optic illumination
- Design allows unhindered access to ear
- Swiveling lens and adjustable speculum holder

Part Number: DH-W35-OL: LED
Part Number: DH-W35-OH: Halogen

Coaxial Ophthalmoscope



- 28 adjustable diopter lenses (+40 to -25)
- 18 unique aperture and filter combinations
- Sealed design keeps dust out
- Anti-theft locking feature
- Halogen and LED available

Part Number: DH-W35-CL-L: LED
Part Number: DH-W35-CH-L: Halogen

Precision Coaxial Ophthalmoscope



- 68 single-diopter steps (+38 to -30) for precise focus
- 18 unique aperture and filter combinations
- Sealed design keeps dust out
- Anti-theft locking feature
- Halogen and LED available

Part Number: DH-W35-PL-L: LED
Part Number: DH-W35-PH-L: Halogen

Apertures



Small, medium, large, slit, fixation

Filters



Cobalt blue, polarization, red-free, unfiltered

Bent Illuminator



- Fiber optic light transmission
- Provides bright consistent illumination

Part Number: DH-W35-BH: Halogen
Part Number: DH-W35-BL: LED

Nasal Illuminator



- Fiber optic light transmission
- Head detaches to convert to throat illuminator

Part Number: DH-W35-NH: Halogen
Part Number: DH-W35-NL: LED

Bi-valve Nasal Illuminator



- Fiber optic light transmission
- Adjustable speculum
- Autoclavable

Part Number: DH-W35-VL: LED
Part Number: DH-W35-VH: Halogen

Throat Illuminator



- Direct light transmission
- Converts diagnostic handle to pen light

Part Number: DH-W35-TL: LED
Part Number: DH-W35-TH: Halogen


Larynx Illuminator



- Fiber optic light transmission
- Rotating wide angle mirror with pilot aperture
- Autoclavable

Part Number: DH-W35-LL: LED
Part Number: DH-W35-LH: Halogen

Diagnostic Station Part Matrix



Diagnostic Accessories

Disposable Ear Specula



- Single use ear specula
- Available in a wide variety of sizes
- Operating otoscope specula also available

DC-SSD-1K-25M: 2.5 mm (0.11"), pack of 1,000
 DC-SSD-1K-40M: 4.0 mm (0.16"), pack of 1,000

Reusable Ear Specula



- Autoclavable up to 134° C (273° F)
- Available in a wide variety of sizes

DC-SSR-01-25M: 2.5 mm (0.11")
 DC-SSR-01-30M: 3.0 mm (0.14")
 DC-SSR-01-40M: 4.0 mm (0.16")
 DC-SSR-01-50M: 5.0 mm (0.20")

Insufflation Bulb



- Otoscope insufflation bulb and tube
- Compatible with fiber optic otoscope

Part Number: DA-OS-IB

Predictive Therm. Probe Covers



- Disposable single-use
- Touchless loading and ejecting of covers

502000: pack of 2000

Infrared Therm. Probe Covers



- Disposable single-use
- Touchless loading and ejecting of covers

Part Number: 303030: Pack of 2,112

Replacement Lamps



- Otoscope & Illuminator
Part Number: DH-W35-BH: Halogen
- Ophthalmoscope
Part Number: DH-W35-BH: Halogen
- Retinoscope
Part Number: DH-W35-BH: Halogen

Two Piece Cuff




- Latex-free bladder
- Variety of sizes available
- Connector options to fit any NIBP device

Stethoscopes



- Stainless steel double-head
- Non-chill rim
- Adult and pediatric available
- Latex-free

Wire Baskets



- 19" x 5" x 3"
(48.26 cm x 12.7 cm x 7.62 cm)
- Various sizes available

Call for Part Number

Distributed By: