

Dual Head Stethoscope



Features:

- Combination diaphragm/Bell chestpiece
- Precision machined stainless steel chestpiece
- Large bell for improved low Frequency response
- Excellent acoustic response of large tunable diaphragm
- Protected rims for patient comfort
- Latex free tubing and eartips
- Provided with extra eartips
- Adult and pediatric options

ST-706BK - Dual Head Adult Stethoscope
 ST-710BK - Dual Head Pediatric Stethoscope

Our dual head stethoscope delivers superior acoustic response from its precision machined stainless steel chestpiece with tunable diaphragm. This design combines superior performance, durability and comfort. This product comes with a 5 year warranty.

5YEAR
warranty



Amico Diagnostic 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada
 Toll Free Tel: 1.877.264.2697 | Tel: 905.763.7778 | Fax: 905.763.8587 www.amico.com