

### Description:

Powerful LED's provide 60,000 Lux in this new innovative exam light. The Gamma 60 LED Light is maintenance free and suitable for any type of exam room. The light comes with four dimming levels, two color temperatures, an antimicrobial finish, and a 360° rotary lighthouse. The Gamma 60 LED light can be ceiling mounted, rail mounted, or on a floor stand to meet all your room requirements.

### Features:

- Rail mount
- Floor Stand
- Ceiling mount
- Wall mount
- Adjustable Color Temperatures
- Adjustable Dimming Levels
- Antimicrobial
- Rotary Feedthrough (360-degree Turnover Axis)



### Specifications:

**Light Head Diameter:** 9.4" (24 cm)

**Light Intensity @ 0.5m:** 60,000 Lux (5,574 fc) at 0.5m (19.7")

**Light Beam Diameter @ 0.5m:** 8.2" (21 cm)

**Color Rendering Index:** 95 Ra

**Color Temperature:** 3500K, 4500K

**Number of LEDs:** 16

**Weight (net):** 9.0 Kg

**Luminaire Body:** Polycarbonate and acrylonitrile butadiene styrene copolymer (PC & ABS)

**Class of Protection:** II

**Electrical:** 100 - 240 VAC, 50/60 Hz, 0.77 amps

**Controls:** Keypad at luminaire head

**Arm System:** Spring balanced articulating arm

**Listings:** Conforms to AAMI STD ES60601-1, IEC STD 60601-1-6, IEC STD 60601-2-41; Certified to CSA STD C22.2 No. 60601-1 by Intertek

**Warranty:** Limited five (5) year warranty

**Ordering Guide:**

**L - GLED60 -**



**For Example: L-GLED60-AB-RH-F**

Gamma Series, Gamma 60 LED Light, Articulating Arm, Rail Mount Horizontal, Fairfield Horizontal Rail

**Ordering Information:**

Series	Arm Type	Mounting Type	Rail Type
L-GLED60	AB Articulating Arm	RH Rail Mount, Horizontal	A No Rail
		RV Rail Mount, Vertical	F Fairfield Horizontal Rail
		RM Rail Mount	X Others
		FX Floor Stand	
		WX Wall Mount	
		CM Ceiling Mount	

**Dimensions:**

