

Description:

The Magna 50 LED magnifier is the lighting solution for the precise examination of parts for a number of fields of application. The Magna 50 LED uses six innovative LEDs, providing critical added value compared to conventional fluorescents magnifiers: four times more light per watt, considerably lowering energy consumption. These LEDs offer a minimum of 50,000 hours of operating life at 100% brightness level – maintenance free for over six years running 24/7.



Specifications:

Connected Load: 100 - 240V; 50/60 Hz

Fitted With: 1 x LED lamp 4W; Neutral White

Work Equipment: Electronic plug-in transformer

Mains Lead: Approx. 7.5 ft (2.3 m)

Luminaire Body Material: Aluminum/plastic

Tubular Sections Material: Steel tube/metal flexible tube

Surface: Painted

Weight (NET): Approx. 2.2 lbs (1 kg)

Fastening:

- Table mount included for table mounted lights
- Rail mount included for rail mounted lights
- Wall mount included for wall mounted lights
- Floor stand included for floor stand light

Technology: Switchable

Usage: Touch-LED

Lamp Cover: Acrylic clear

Class of Protection: II

Color Luminaire Body: Black

Color of Tubular Sections: Silver metallic

Lens Dimensions: Ø 2.8" (72 mm)

Lens Diopters: 5.3

Add-On Lens: Not available

Ordering Guide:

L - MA50 -



Model Examples:

L-MA50-GB-RH-A	Magna 50, Gooseneck Arm, Basic, Rail Mount, Horizontal, Amico Rail
L-MA50-GB-WX	Magna 50, Gooseneck Arm, Basic, Wall Mount
L-MA50-GE-FX	Magna 50, Gooseneck Arm with Extension, Floor Stand
L-MA50-GE-TX	Magna 50, Gooseneck Arm with Extension, Table Mount

Ordering Information:

Series	Arm Options	Mounting Options	Rail
L-MA50	GB Gooseneck Arm, Basic	RH Rail Mount, Horizontal	A Amico Rail (For Rail Mount Only)
	GE Gooseneck Arm with Extension	RV Rail Mount, Vertical	F Fairfield Rail (For Rail Mount Only)
		FX Floor Stand	X Not Applicable
		WX Wall Mount	
		TX Table Mount	

***NOTE:**

North American: 120 V; International: 230 V

Dimensions:

