

# Magna Series

## Magna 9 Magnifier

### Description:

Magna 9 features a distinctive fully rotational 3-diopter magnifier lens for precise positioning. Plus three 9-Watt fluorescent bulbs provide two settings for adjusting the light output to 18 or 27 watts.



### Specifications:

**Connected Load:** 120 V; 60Hz

**Fitted With:** Approx. 2 m cordset

**Work Equipment:** Conventional ballast

**Mains Lead:** Approx. 9.8 ft (3 m); shock-proof plug

**Luminaire Body Material:** Plastic

**Tubular Sections Material:** Steel Tube

**Surface:** Painted

**Weight (NET):** Approx. 9.9 lbs (4.5 kg)

#### Fastening:

- Table mount included for table mounted lights
- Rail mount included for rail mounted lights
- Wall mount included for wall mounted lights
- Floor stand included for floor stand light

**Technology:** Multi-Stage Switchable

**Usage:** Switch

**Lamp Cover:** Screen; Milky

**Class of Protection:** I

**Color Luminaire Body:** Light grey

**Color of Tubular Sections:** Light grey

**Lens Dimensions:** 162 x 105

**Lens Diopters:** 3

**Special Features:** Light output controlled with 2-way switch

# Magna Series

## Magna 9 Magnifier

### Ordering Guide:

**L - MA09 -**



### Model Examples:

<b>L-MA09-AB-RV-F</b>	Magna 9, Articulating Arm, Basic, Rail Mount, Vertical, Fairfield Rail
<b>L-MA09-AB-WX</b>	Magna 9, Articulating Arm, Basic, Wall Mount
<b>L-MA09-AE-FX</b>	Magna 9, Articulating Arm with Extension, Floor Stand
<b>L-MA09-AE-TX</b>	Magna 9, Articulating Arm with Extension, Table Mount

### Ordering Information:

Series	Arm Options	Mounting Options	Rail
<b>L-MA09</b>	<b>AB</b> Articulating Arm, Basic	<b>RH</b> Rail Mount, Horizontal	<b>A</b> Amico Rail (For Rail Mount Only)
	<b>AE</b> Articulating Arm with Extension	<b>RV</b> Rail Mount, Vertical	<b>F</b> Fairfield Rail (For Rail Mount Only)
		<b>FX</b> Floor Stand	<b>X</b> Not Applicable
		<b>WX</b> Wall Mount	
		<b>TX</b> Table Mount	

### \*NOTE:

North American: 120 V; International: 230 V

### Dimensions:

