

GoLift

Specifications

Product Features:

- Wired Hand Control
- Quick Connect Carry Bar
- Quick Connect Hand Control
- Light Weight and Portable
- LED Indicators
- Infection Control
- Durable and Strong
- Easy to Service

Safety Features:

- Emergency Stop
- Emergency Lowering Device
- Upper Limit Detection
- Lower Limit Detection
- Slack Tape Sensor
- Free Fall Brake (over speed governor)
- Low Battery and Dead Battery Alarms
- Soft Start and Stop
- Overload Protection
- Emergency Manual Lowering



Single Lift GO-LFT-400

GoLift 400
Lifts up to 400 lbs / 181 kg



Single Lift GO-LFT-700

GoLift 700
Lifts up to 700 lbs / 318 kg



Tandem Lift GO-LFT-1000

GoLift 1000
Lifts up to 1000 lbs / 454 kg



Specifications:

Maximum Lifting Speed

No load: 2 inches/second
 GoLift400: 1 inch/second
 GoLift700: 1 inch/second
 GoLift1000: 1 inch/second

Unit Weight

GoLift400: 10 lbs (4.5 kg)
 GoLift700: 10 lbs (4.5 kg)
 GoLift1000: 20 lbs (9 kg) *(Tandem Lifts)*

Safe Working Load (SWL)

GoLift400: 400 lbs (181 kg)
 GoLift700: 700 lbs (318 kg)
 GoLift1000: 1000 lbs (454 kg)

Dimensions

GoLift400 / GoLift700:
 • Length: 7.75" (196.8 mm)
 • Width: 7.75" (196.8 mm)
 • Height: 4.25" (107.9 mm)
 • Strap Length: 84" (2133.6 mm)

GoLift1000:
 • Length: 15.5" (393.7 mm)
 • Width: 7.75" (196.8 mm)
 • Height: 4.25" (107.9 mm)
 • Strap Length: 84" (2133.6 mm)

Number of Lifts Per Charge (Duty: 10/90)

25% of Strap at Midrange¹:
 • 375 with 185 lbs (84 kg)
 • 170 with 400 lbs (181 kg)
 • 102 with 700 lbs (318 kg)

Charging time: 2 - 4 hours
¹Calculated using 5Ah Packs

Lift Case

Flame retardant ABS

Batteries

High Capacity, Nickel Metal Hydride (NiMH)
 Standard: 2x 14.4V (3.3Ah)
 Optional: 2x 14.4V (5Ah)

Battery Charger

Input: 100-240V AC
 Output: 36V DC, 1.0A, 40W

Hand Control

Protection Class: IPX4

Maximum Sound Level

<55 dB


Expected Service Life

10 years or 22,500 cycles


Approvals

Certified to: CAN/CSA-C22.2
 No 60601-1:08
 UL 60601-1:08
 Tested to: ISO 10535-06, CE


GoTrack Configurations:



Regular



Double



Triple

