

# Patient Lift Slings



Amico



# Amico Patient Lift Slings

Amico Mobility Solutions is dedicated to manufacturing products with quality and comfort in mind for both patients and clinicians. Amico's patient slings are available in various sizes, colors, materials, and customizations to meet every need, including colored sling loops for ease of identification.

## Our Slings

**Go Basic**

**Go Hygiene**

**Go Comfort**

**Go Amputee**

**Go Long**

**Go Active**

**Go Reposition**

**Go Turn**

**Go Walk**

**Go Limb**

Visit Amico's website [www.amico.com/videos](http://www.amico.com/videos) to see our slings in action!











## Application

- Full body with optional head support
- Lift and transfer patient from any surface including the floor
- Reusable slings lift up to 1,000 lbs (454 kg)
- Disposable slings lift up to 660 lbs (300 kg)

## Available Materials

● Reusable    ● Disposable

## Best Used For



Patient transfer



Boosting







## Application

- Assistance to clients who require hip and thigh support
- Support for patients with amputated limbs
- Lift and transfer or position patients from seated or laying position
- Support for patients with head control but reduced upper body strength
- Reusable slings lift up to 1,000 lbs (454 kg)
- Disposable slings lift up to 660 lbs (300 kg)
- GoComfort can also be used as an amputee sling



## Available Materials

● Reusable
 ● Disposable

## Best Used For



Amputee support





## Application

- Support for the transfer of a patient to other surfaces including beds and stretchers
- Assistance to reposition the patient on bed
- Reusable slings lift up to 1,000 lbs (454 kg)
- Disposable slings lift up to 660 lbs (300 kg)

## Available Materials

● Reusable    ● Disposable

## Best Used For



Supine/prone positioning



Supine/lateral positioning



Supine/supine transfer





# Go Reposition



## Application

- Assistance to reposition the patient on a bed
- Additional positioning slings included
- Reusable slings lift up to 300 lbs (136kg)

## Available Materials

● Reusable

## Best Used For



Prone/lateral positioning



Prone/supine positioning







## Application

- Walking and balance training support
- Designed for patients who need support for ambulation
- Additional trunk and hip support from the Amico GoActive sling
- Disposable slings lift up to 660 lbs (300 kg)

## Available Materials

● Reusable    ● Disposable

## Best Used For



Support for gait  
training and  
ambulation







## Application

- Assistance with toileting and hygiene through large sling opening
- Support for patients with head control but reduced upper body strength
- Adjustable belt for support around the body
- Reusable slings lift up to 1,000 lbs (454 kg)
- Disposable slings lift up to 660 lbs (300 kg)

## Available Materials

● Reusable    ● Disposable

## Best Used For



Toileting



Bathing



Peri-care







## Application

- Support for lifting bilateral and unilateral amputees
- Disposable slings lift up to 660 lbs (300 kg)


## Available Materials

● Disposable

## Best Used For



Amputee support







## Application

- Walking and balance training support
- Special support for patients with reduced movement abilities
- Reusable slings lift up to 600 lbs (272 kg)
- Disposable slings lift up to 660 lbs (300 kg)

## Available Materials

● Reusable    ● Disposable

## Best Used For



Support for gait  
training and  
ambulation







## Application

- Provides sling solutions for repositioning
- Reusable slings lift up to 600 lbs (272 kg)
- Disposable slings lift up to 660 lbs (300 kg)

## Available Materials

● Reusable    ● Disposable

## Best Used For



Prone/lateral positioning







## Application

- Support for lifting patient arms or legs
- Support for passive range of motion exercises
- Reusable slings lift up to 600 lbs (272 kg)
- Disposable slings lift up to 660 lbs (300 kg)

## Available Materials

- Reusable
 ● Disposable

## Best Used For




Support for limbs
















# Sling Details and Sizing

## Guidelines For Choosing Your Sling Sizes



Slings	S	M	L	XL
 With Head Support	Height: 29" Shoulder: 20" Hip: 16"	Height: 30" Shoulder: 24" Hip: 20"	Height: 31" Shoulder: 28" Hip: 24"	Height: 32" Shoulder: 32" Hip: 28"
 No Head Support	Height: 21" Shoulder: 20" Hip: 16"	Height: 22" Shoulder: 24" Hip: 20"	Height: 23" Shoulder: 28" Hip: 24"	Height: 24" Shoulder: 32" Hip: 28"
	Height: 23" Shoulder: 24" Hip: 20"	Height: 24" Shoulder: 26" Hip: 22"	Height: 25" Shoulder: 28" Hip: 24"	Height: 26" Shoulder: 30" Hip: 26"
	42" x 78"		56" x 78"	
	20" x 75"			
	Hip: 15"	Hip: 18"	Hip: 22"	Hip: 26"
	Chest: 32"	Chest: 36"	Chest: 40"	Chest: 44"
	Height: 28" Hip: 15"	Height: 33" Hip: 18"	Height: 38" Hip: 22"	Height: 43" Hip: 26"
	Chest: 40"	Chest: 44"	Chest: 48"	Chest: 52"
	Height: 31" Width: 18.5"	Height: 35" Width: 19.5"	Height: 39.5" Width: 19.5"	Height: 42.5" Width: 21"
	15" x 7"			

## Material Information


	Mesh	Polyester	Poly Woven
Reusable	◆		◆
Disposable		◆	



\*Includes manufacturing date\*



# Sling Best Use Cases



	Basic	Comfort	Long	Reposition	Walk	Hygiene	Amputee	Active	Turn	Limb	Custom
Patient transfer	♦	♦	♦			♦	♦				♦
Boosting	♦		♦								♦
Supine/prone positioning			♦	♦					♦		♦
Supine/lateral positioning			♦	♦					♦		♦
Supine/supine transfer			♦								♦
Toileting, bathing and peri-care						♦					♦
Amputee support		♦					♦				♦
Support for gait training					♦			♦			♦
Support for limbs										♦	♦
Prone/lateral positioning				♦					♦		♦

## Washing Instructions for All Slings

Normal wash at indicated temperature	Do not bleach	Tumble dry low	Do not iron	Do not dry clean





Amico Mobility Solutions Corporation | 122-B East Beaver Creek Road, Richmond Hill, ON L4B 1G6, Canada  
600 Prime Place, Hauppauge, NY 11788, USA  
Toll Free Tel: 1.877.462.6426 | Tel: 905.764.0800 | Fax: 905.764.0862  
Email: [info@amico.com](mailto:info@amico.com) | [www.amico.com](http://www.amico.com)