

Amico 4 Port Turntable Specifications



General Specifications:


The Amico 4 Port Turntable is suitable for installations where junctions are required in a track or when a track turn is required that a standard curved sections cannot achieve. The Turntable provides ultimate flexibility in lift and transfer situations that demand a wide range of locations, or require specialized options.

The four (4) exits of the Turntable allow multiple configuration options for the system designer and gives the user more options for the installation. An adaptation that is being planned can be designed around the angles of the turntable. For example, a turntable in a bathroom could be positioned so as to allow transfer between chair, bath, wheelchair, shower area and or changing stretcher.

Diagrams



Typical Installation using Unistrut hardware and 1/2" (12.7 mm) threaded rods



Location of 4 x 1/2" (114.3 mm) holes for the 1/2" (12.7 mm) threaded rods