

Baby Bassinet

Marco Series



Configuration Examples








Customize Your Bassinet



Options

				
Pull Out Writing Surface	Flip Up Writing Surface	Flame Retardant Mattress (26" L x 12" W x 2" H)	Infant Tub (26.5" L x 13" W x 10" H)	Drawer Dividers
				
Laminate Chart Holder	Plastic Chart Holder	Key Lock	Punch Code Lock	Electronic Lock

Standard Finishes*

						
Hard Rock Maple	Candlelight	Dolce Vita	Medium Oak	Chocolate Pear	Snow White	Richmond Grey

*Custom colors are available, please contact us for details.

Available Materials

Melamine, Laminate, or Thermofoil

Standard Dimensions**

Overall Height (floor to handle)	37" (94 cm)
Height from Floor to Top of Mattress (floor to baby)	34.25" (87 cm)
Overall Width	32" (81.3 cm)
Overall Depth	21" (53.3 cm)
Casters	5" (12.7 cm) 2 locking, 2 swivel
Flip Up	14" L x 17" D (35.6 cm L x 43.2 cm D)

**Custom dimensions are available, please contact us for details.