

Bedside Cabinets



BC-JD3-0BB-1CCL



BC-JDD-0XB-1CHM

Jasper Series

- High pressure laminate top and drawer fronts
- Casters (exposed or recessed) or self-levelling glides
- Different drawer configurations
- Seamless spillguard top available

Oliver Series

- Thermofoil top & drawer fronts
- Seamless for infection control
- Casters (exposed or recessed) or self-levelling glides
- Different drawer configurations
- Seamless spillguard top available



BC-ODD-1BT-2GCP



BC-OD3-1BA-4CHM



BC-LDD-0B-C1MC-L



BC-LDO-0B-G1HM-L

Lucas Series

- Anodized aluminum legs for durability
- Laminate or thermofoil tops and drawer fronts
- Casters (exposed) or self-levelling glides
- Different drawer configurations
- Seamless spillguard top available

Customize Your Bedside Cabinets



Drawer Configurations



3 DRAWERS



DRAWER/DOOR



DRAWER/OPEN

Features and Options



Steel Side Drawer Construction




Vented Back



Casters
(2 locking, 2 unlocking)



Seamless Spillguard Top



Field Replaceable Drawer and Door Fronts

Drawer/Door Fronts



Flat



Bevelled*



Convex*



Shaker*

*Only available in Thermofoil material.

Hardware



Modern Stainless Steel Pull



Stainless Steel Knob



Brushed Nickel Pull



Modern Brass Pull



Brass Knob


Standard Finishes**



Hard Rock Maple




Candlelight



Dolce Vita



Medium Oak



Chocolate Pear



Snow White



Richmond Grey

**Custom colors are available, please contact us for details.

Standard Dimensions

	Height	Width	Depth
Jasper Series	32"	20"	21"
Oliver Series	32"	20"	21"
Lucas Series	32"	20"	21"