


Overbed Tables



Various top styles available



Spillguard top holds up to 36 oz

Vanity drawer with mirror



Seamless cup holders




Easy, one-hand activated secondary surface




Integrated accessory rail (urinal hook also available)

Shock-absorbent gas spring cylinder for easy height adjustment



OT-F-TYYS-XU

Various base styles available



Lowest base on the market
Exposed steel; no plastic cover allows for easy cleaning

Configuration Examples

Single Top



OT-S-TYVV-XU

Dual Top



OT-F-TYYS-XU

Specialty



JELLYBEAN SHAPE



TABLET MOUNT



OT-S-TYVN-XC



OT-F-TYYS-XU



WHITEBOARD



SOLID SURFACE

Customize Your Overbed Tables



Table Top Options

Standard (Laminate)	Standard (Thermofoil)	Specialty (Thermofoil only)	Corian/Solid Surface	Kidney Shape (32" x 19")	Peanut Shape (32" x 19")	Jellybean Shape (35" x 22")	Whiteboard (Laminate only)

Standard Top Finishes*

Hard Rock Maple	Candlelight	Dolce Vita	Medium Oak	Chocolate Pear	Snow White	Richmond Grey

Standard Base Finishes*

Graphite Grey	Nevada	Almond

*Custom finishes are available, please contact us for details.

Accessories

Vanity Drawer and Mirror	Urinal Hook	Locking Casters	Storage Baskets (10" x 6" x 9" and 5" x 5" x 4")

Flip Table Tops**

Flat Surface	Dual Cup Holders	Spillguard

**Flips are not available on laminate, storage, jellybean, c-base, or solid surface overbed tables.

Storage Table Options

Shelf 1 (with 2 compartments) Overall Dimensions: 17" L x 10" W x 1.5" D	Shelf 2 (with 1 compartment) Overall Dimensions: 17" L x 10" W x 1.5" D	Storage Drawer (with 1 compartment) Overall Dimensions: 17" L x 10" W x 4" D

Storage Table Configuration Examples

OT-ST-TYYN-XU-E	OT-ST-TYYV-XU-E

Standard Dimensions**

Single Table Top	32" L x 17" W (81.3 cm x 43.2 cm)
Flip Table Top	13.5" L x 17" W (34.3 cm x 43.2 cm)
Storage Table Top	43.5" L x 19" W (110.5 cm x 48.3 cm)
Height Range	28" Low - 42" High (71 cm - 106.7 cm)

**Custom dimensions are available, please contact us for details.

	U Base	H Base	C Base
Overall Length	32.5"	34.38"	30.75"
Overall Width	19"	18.25"	21.25"
Internal Width (between the legs)	14"	29.38"	25"
Overall Height from the Floor	Approx. 1.9"	Approx. 1.9"	Approx. 1.9"