

Tablet Mount Overbed Tables



Tablet holders available for:

- iPad
- iPad Mini
- Windows Surface
- Samsung Galaxy



Tablet locks



Flip surface option




USB ports

Power cord



Two handles for height adjustment

Integrated accessory rail
(urinal hook and small basket)




Cord wrap

TM-S-TYYN-GU-A

14"

Various base styles available



U Base








H Base

Lowest base on the market
Exposed steel; no plastic cover allows for easy cleaning


Customize Your Overbed Tables






Table Top Options

| | | | | | | |
|---|---|---|---|--|---|---|
|  |  |  |  |  |  |  |
| Standard (Laminate) | Standard (Thermofoil) | Specialty (Thermofoil only) | Corian/Solid Surface | Kidney Shape (32" x 19") | Peanut Shape (32" x 19") | Whiteboard (Laminate only) |

Standard Top Finishes*

| | | | | | | |
|---|---|---|---|---|---|--|
|  |  |  |  |  |  |  |
| Hard Rock Maple | Candlelight | Dolce Vita | Medium Oak | Chocolate Pear | Snow White | Richmond Grey |

Standard Base Finishes*



| | | |
|---|---|---|
|  |  |  |
| Graphite Grey | Nevada | Almond |

*Custom colors are available, please contact us for details.

Optional Accessories


| | | | |
|--|--|--|--|
|  |  |  |  |
| Vanity Drawer and Mirror | Urinal Hook | Locking Casters | Storage Baskets (10" x 6" x 9" and 5" x 5" x 4") |

Flip Table Tops**



| | | |
|---|--|--|
|  |  |  |
| Flat Surface | Dual Cup Holders | Spillguard |

**Flips are not available on laminate, or solid surface overbed tables.

Tablet Mounting Option

| |
|---|
|  |
| Arm |

USB Ports Options

| | |
|---|--|
|  |  |
| 4 x USB-A | 2 x USB-A + 1 USB-C |

Standard Dimensions**

| | |
|------------------|---------------------------------------|
| Single Table Top | 32" L x 17" W (81.3 cm x 43.2 cm) |
| Flip Table Top | 13.5" L x 17" W (34.3 cm x 43.2 cm) |
| Height Range | 28" Low - 42" High (71 cm - 106.7 cm) |

**Custom dimensions are available, please contact us for details.

| | U Base | H Base |
|-----------------------------------|--------------|--------------|
| Overall Length | 32.5" | 34.38" |
| Overall Width | 19" | 18.25" |
| Internal Width (between the legs) | 14" | 29.38" |
| Overall Height from the Floor | Approx. 1.9" | Approx. 1.9" |