

Recliner Chairs



Scarlett Series:

- 3 position recliner chair (24°, 35°, and 52° recline)
- Available in standard and bariatric models
- Ideal for visitors and short term stay patients
- Options include: IV pole, arm caps, flip up tables, and heat and massage

Daphne Series:

- Infinite recline positions from seated to fully flat
- Available in standard and bariatric models
- Integrated recessed handle to activate recline
- Pull-lever to activate foot rest
- Extra wide seat and arms ideal for LDR and NICU areas
- Options include: Trendelenburg, IV pole, arm caps, foot transport tray, flip up tables, and heat and massage



Customize Your Recliner Chairs



Features and Options - Scarlett Series

				
Polyurethane or Wood Arm Caps	IV Pole	3" Non-Marking Casters	Flip Up Side Table	Pivoting Table (Single or Double)
				
Metal Push Bar	Pillow Back	Central Braking	Heat and Massage	Accessory Holder

Features and Options - Daphne Series

						
Trendelenburg	IV Pole	Lever-Controlled Foot Rest	Recessed Handle to Recline	Polyurethane or Wood Arm Caps	Metal Push Bar	Pillow Back
						
Central Braking	3" Non-Marking Casters	Pivoting Table (Single or Double)	Flip Up Side Table	Pull Out Foot Tray	Heat and Massage	Accessory Holder

Standard Dimensions

Description	Scarlett Series Standard	Scarlett Series Bariatric	Daphne Series Standard	Daphne Series Bariatric
Overall Dimensions	46.75" H x 29.25" W x 37.5" D	46.75" H x 33" W x 37.5" D	44" H x 33" W x 39" D	44" H x 37" W x 39" D
Seat Dimension	20" W x 21" D x 19" H	24.5" W x 21" D x 19" H	23.5" W x 21" D	28.5" W x 21" D
Depth Reclined	70"	70"	73"	73"
Weight Capacity	350 lbs	500 lbs	500 lbs	750 lbs