





LAB ROTARY VANE - DRY DUPLEX TANK MOUNTED VACUUM SYSTEM (0.5 HP - 8.9 HP)



DIMENSIONS SUBJECT TO CHANGE WITHOUT NOTICE

*SYSTEM CONFIGURATION DIFFERS FROM MODEL SHOWN

**380 VOLT DATA PROVIDED AT 50HZ
ALL OTHER VOLTAGES ARE RATED AT 60HZ

ADDITIONAL MODELS AVAILABLE UPON REQUEST

MODEL	HP (KW)	A INLET NPT	B OUTLET NPT	TANK SIZE GALLON (LITRE)	PUMP SCFM @ 19" HG (LPM)		SYSTEM CAPACITIES @ 19" HG (LPM)		WEIGHT IN LBS (KG)	SQ. FT. REQUIRED (M ²)	SYSTEM FLA (WITH ALL PUMPS RUNNING)				SYSTEM BTU/HR	SOUND LEVEL dB (A)	W (M)	L (M)	H (M)
					50 Hz Motor	60 Hz Motor	50 Hz Motor	60 Hz Motor			208V	230V	380V**	460V					
V-RVD-D-080P-TH-L-005*	0.5 (0.37)	0.5	1.0	80 (303)	2.0 (57)	2.4 (68)	4.0 (113)	4.8 (136)	685 (311)	18.8 (1.75)	6	6	4	3	2,544	65	36 (0.91)	75 (1.91)	60 (1.52)
V-RVD-D-080P-TH-L-007*	0.75 (0.56)	0.5	1.0	80 (303)	2.8 (79)	3.4 (96)	5.6 (159)	6.8 (193)	703 (319)	18.8 (1.75)	9	8	5	4	3,816	66	36 (0.91)	75 (1.91)	60 (1.52)
V-RVD-D-080P-TH-L-012	1.2 (0.89)	0.75	1.0	80 (303)	4.4 (125)	5.3 (150)	8.8 (249)	10.6 (300)	727 (330)	19.8 (1.84)	12	12	7	7	6,106	69	38 (0.97)	75 (1.91)	60 (1.52)
V-RVD-D-080P-TH-L-017*	1.7 (1.27)	0.75	1.0	80 (303)	6.7 (190)	8.0 (227)	13.4 (379)	16.0 (453)	781 (354)	19.8 (1.84)	16	15	8	8	8,652	73	38 (0.97)	75 (1.91)	60 (1.52)
V-RVD-D-080P-TH-L-032*	3.2 (2.39)	1.25	1.25	80 (303)	11.3 (320)	13.5 (382)	22.6 (640)	27.0 (765)	961 (436)	20.8 (1.95)	20	18	9	9	16,284	79	40 (1.02)	75 (1.91)	60 (1.52)
V-RVD-D-080P-TH-L-046*	4.6 (3.43)	1.25	1.25	80 (303)	14.2 (402)	17.0 (481)	28.4 (804)	34.0 (963)	1005 (456)	22.4 (2.08)	-	27	13	13	23,408	80	43 (1.09)	75 (1.91)	64 (1.63)
V-RVD-D-120P-TH-L-056*	5.6 (4.18)	1.25	1.25	120 (454)	17.5 (496)	21.0 (595)	35.0 (991)	42.0 (1189)	1272 (577)	30.6 (2.84)	33	32	17	16	28,498	79	55 (1.40)	80 (2.03)	64 (1.63)
V-RVD-D-120P-TH-L-089*	8.9 (6.64)	1.5	1.5	120 (454)	23.3 (660)	28.0 (793)	46.6 (1320)	56.0 (1586)	1340 (608)	30.6 (2.84)	51	49	25	25	45,290	80	55 (1.40)	80 (2.03)	64 (1.63)



85 Fulton Way, Richmond Hill
Ontario, L4B 2N4 Canada
71 East Industry Court
Deer Park, NY 11729
Direct Tel: 905.764.0800
Tel: 1.877.462.6426
Fax: 1.866.440.4986